
# ACUAVEST®

HOJA DE INSTRUCCIONES	ES	2
WORKING INSTRUCTIONS	EN	4
ISTRUZIONI PER L'USO	IT	6
NOTICE D'EMPLOI	FR	8
INFORMACIÓN ADICIONAL ADDITIONAL INFORMATION INFORMAZIONI SUPPLEMENTARI INFORMATIONS COMPLÉMENTAIRES		10

Consulte nuestra web  
Visit our website

# ACUAVEST®

## DESCRIPCIÓN

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# ACUAVEST®

## DESCRIPTION

Special investment for Cobalt Chrome alloys. According to ISO 15912:2016, the classification of the product is Type 2 Class 1.

This material contains quartz. **DO NOT INHALE THE POWDER!**

## INSTRUCTIONS OF USE

ACUAVEST® is designed to be used with running water. However, a mixing liquid **EXPANSOR** must be used in the preparation of the model.

**DO NOT USE EXPANSOR LIQUID WHICH HAS BEEN FROZEN (BECOME CRYSTALLISED).**

### MIX PROPORTION

#### For the duplicate model:

Use 200 g of investment powder and mix with the following liquid:

13-14 ml EXPANSOR + 13-14 ml distilled water

#### For the cylinder:

Use 400 g of investment powder and mix with the following liquid:

55-60 ml distilled water

In both cases, if the mixture is made by hand, mix it with a vibrator. Put the liquid in first, adding the powder afterwards and beating for 60 seconds. If mixing under vacuum, mix for 30 sec.

**WARNING: DO NOT USE CONTAINERS THAT HAVE BEEN IN CONTACT WITH PLASTER.**

### WORKING TIME

#### **4 or 5 minutes.**

Room temperature, powder & liquid at 20-22°C.

If the liquid or gelatines are colder the setting time is longer.

### DUPLICATION OF MODELS

The duplication can be carried out in gelatin (TECHNOJEL®) or silicone (DUPLIFLEX).

#### FILLING THE MODELS

Pour the mixed investment into the duplicate under vibration. To guarantee a smooth and hard surface of the models **it's important to maintain the vibration for 2 minutes.**

#### SETTING TIME

**30-35 minutes**, both for models as well as cylinders.

#### DRYING THE MODELS

**40 minutes** in a furnace at **220-240°C.**

#### MODEL HARDENING


After drying the models they may be hardened using any system. We recommend our hardener TECHNODUR. If silicon moulds are used, no hardening is necessary.

## PREPARATION FOR INVESTING

- Preferably use fine investments for the first coat (TECHNOPAINT) or degrease the wax sufficiently (TENSIO-WAX).
- Fill in the cylinder with a low vibration and **demould at 35 minutes**.

## HEATING

### DRYING

After setting for 35 minutes put the cylinders in the furnace, preferably cold, and heat to 240 °C for 40 minutes.


### HEATING

Once the cylinders are dry, programme the casting temperature.

We recommend no more than 900°C for partial pieces and 1.000°C for whole plates.

Maintain this temperature for at least 35 minutes.

The best expansion rates are obtained at a heating of 5°C per minute.

If the indicated drying times and pre-heating temperatures are not observed, slight fissures in the muffle surface may occur; they are of no importance for the cast object.

## AFTER CASTING

Wait until the cast metal gets to black.

Cool the cylinders with water. Act carefully: firstly, wet slowly and when the cylinder is colder, then submerge it completely in water. This makes easier the removing of the moulds.

## STORAGE AND EXPIRY


Batch number and expiry date are marked on each package.

### POWDER ACUAVEST®

Store at 18-25°C.

Keep the packages tightly closed in a dry place.

The investment is sensitive to moisture: prolonged exposure to the air may alter the characteristics of the product.

### EXPANSOR

Store at 18-25°C.

Prolonged exposure to temperatures lower than 5°C or higher than 35°C may produce freezing and/or destabilisation of the product.

Do not use the liquid if it shows a sign of cristallization.

Keep the bottles tightly closed after each use. Prolonged exposure to the air may cause crystallisation.

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EN

Our recommendations referring to the techniques applied, whether verbally, in writing or in the form of practical instructions, are based on our own experiences and/or tests, and must therefore be considered only as instructions.

For further information, see the safety data sheet of ACUAVEST® and EXPANSOR. Order them to your local dealer or visit our website to download them.

IMPORTANT: Check our website for the most updated Instruction Manual version, as this information is subject to continuous review.

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