R-100118

E.....	CÁMARA DE MICROCHORREADO RED.....	3
GB.....	ROUND MICRO-BLASTING CHAMBER.....	4
D.....	RUNDE FEINSTRAHLKAMMER.....	5
P.....	CÂMARA DE MICRO JACTO REDONDA.....	6
I.....	CAMERA DI MICROSABBIATURA TONDA.....	7

MESTRA®


TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

T (+34) 944530388 - F (+34) 944711725

E-mail: mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

- | | |
|---|---|
| A | Ventana / Window / <i>Fenster</i> / Janela / Finestra |
| B | Aberturas laterales / Side openings / <i>Seitliche Öffnungen</i> / Aberturas laterais / Aperture laterali |
| C | Interruptor / Main switch / <i>Schalter</i> / Interruptor / Interruttore |
| D | Fusible / Fuse / <i>Sicherung</i> / Fusível / Fusibile |
| E | Conmutador tensión / Voltage switch / <i>Spannungsumschalter</i> / Comutador de tensão / Commutatore tensione |
| F | Ventilador / Fan / <i>Lüfter</i> / Ventilador / Ventilatore |
| G | Bandeja para polvo / Dust tray / <i>Staubwanne</i> / Bandeja para pó / Raccogliore polvere |
| H | Lengüeta / Tab / <i>Zunge</i> / Lingueta / Linguetta |

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The R-100118 micro-blasting chamber is a compact and light dust collector, optimal for a small dental lab. It is a perfect partner for the R-100122 Micro-blast Unit (a practical mini-blasting unit specially designed for the dental sector), and for collecting abrasives and sands after the blasting, cutting and drilling processes.

The unit features a silent and efficient extractor fan system. Its wide transparent window can be dismantled. The abrasive is collected in a lower tray (G) and therefore it can be easily removed.

INSTALLATION

Remove the micro-blasting chamber from its packing. The transparent lid (A) can be separated from the base by pulling it upwards. To place it again it must be centred with a small protrusion in the base.

Select the suitable voltage for your power supply (220/127 V) by pressing the voltage switch (E).

220 V

127 V

Start the unit by turning the main switch (C) on. The fan (F) will start up. While spinning the fan will create a negative pressure inside the chamber, and the abrasives and sands will be kept inside the unit. Make sure that no object or material interferes in the free air flow from the fan; otherwise the unit could overheat or malfunction.

OPERATION

The working piece is inserted into the chamber through one of the side openings (B), while the micro-blasting unit or drill motor is inserted through the opposite side. The piece or equipment does not always need to be inserted completely in the unit: it is enough to insert the part of the piece or equipment needed to carry out the work. This is learnt through experience and comfort. The progress of the work can be observed through the transparent window (A).

DO NOT AIM THE BLAST UNITS STRAIGHT AT THE FAN.

Upon finishing the work, wait a few seconds for the abrasives and dust to settle before taking your hands or the pieces out of the chamber.

Always stop the unit after using it. It will prevent the fan working unnecessarily.

To remove the dust collected in the lower tray, pull the tab (H) and carefully open the base.

MAINTENANCE AND CLEANING

To ensure a long life of the unit and an optimal operation, we strongly recommend:

- To use the unit only in networks with specially treated air. Humidity and oil should be specially avoided.
- To use a high quality abrasives. An inappropriate grain will cause the unit wearing out prematurely in the medium-term. The used abrasive should not be recycled.
- To make always sure that the unit has been turned off and unplugged before dismantling or cleaning the chamber.
- To clean the chamber, wipe it with a cloth.
- To wipe the transparent lid (A) and its inner surface smoothly with a damp cloth to remove the abrasives or sands. Then place them again in the base.

CAUTION

- Do not ever exceed 6 bar pressure . The optimal working pressure lies between 4 and 5 bars.
- Work correctly protected: always wear gloves, mask and protection goggles.
- Do not allow children and non-qualified personnel to handle the unit.
- Only use high quality abrasives (it is strongly recommended the use of **MESTRA** abrasives). Poor quality abrasives could be harmful by inhalation and can obstruct the nozzles.
- Make sure the pressure air installation to be fitted to the unit is correctly protected by at least one pneumatic discharge valve.
- In case of leaks along the pressure air circuit, consult the installer immediately.

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