## R-100136

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R-100136 is a practical Ultrasonic Bath specially designed for a wide range of applications. For a long working life and an optimum performance, we recommend you to read carefully the following instructions:

## INSTALLATION

- Carefully unpack the Ultrasonic Bath. Inside you will find the following accessories:
  - 1 draining hose (M).
  - 1 power cord (L).
  - 1 plastic basket (K).
- Choose a flat and firm surface far from heat or vibration sources.
- Connect the hose (M) to the drainage pipe (D).
- Connect the power cord (L) to the socket of the unit (J) and to a 230 V, 50 Hz earthed power source.
- Open the lid (A), place the products in the basket (K) in line and then put the basket inside the tank (B).
- In order to prevent scratches and friction of the products, do not pile them up.
- Pour water into the stainless steel tank. The water should cover the products but within the *min.* and *max.* marks.
- The machine is ready for operation.

## CONTROL PANEL AND OPERATIONS

### 1. Working time display


It counts down after work is started.

### 2. Multistage thermometer:

Green	Water temperature is between 0 and 30 % of set temperature.
Orange	Water temperature is between 30 and 60 % of set temperature.
Red	Water temperature is between 60 and 90 % of set temperature.
Flashing red	Water temperature is close to the set temperature.
All lights	Water temperature is at the set temperature (90-100 %)

### 3. Set temperature display

5 set temperatures can be selected. Press  button to select the temperature in the following sequence:

### 4. °C or °F indicator

### 5. Solution usage remaining time

It counts down to show the remaining usage time (hour:minute) of the solution.

### 6. Solution warning time

Once the solution usage timer is set, the light bulb icon comes up.

 flashing, indicates that solution needs to be changed.

### 7. Degas

Illuminated: The degas function is selected.

Flashing: Degassing is underway.

### 8. Normal cleaning status

Illuminated: Normal cleaning is underway.

### 9. Heater status

Press  button **twice**. The icon  flashes, indicating the heater is turned on. Press  button **twice** again, the icon  dims, indicating the heater is turned off.

## 10. Working status

Green: Normal working status, ready to use.

Red: Under overheat protection. 20 minutes or more of waiting time is needed for the light to change from red to green. Work can resume afterwards.

## 11. Heater button

Press  button **twice**, the icon  flashes, indicating the heater is turned on. Press  button **twice** again, the icon  dims, indicating the heater is turned off.

When the water reaches the set temperature, the heater will turn off automatically.

## 12. Function button

Press **Function** button . The icon  illuminates, indicating the heater is turned on. To cancel heating, press  button again,  dims, indicating the heater is turned off.

### Degas:

When cleaning heavily soiled or greasy items or silver, copper or brass items, chemical solutions need to be added into the water. The solution may form many small droplets and take a long time to dissolve in the water. Newly added water may generate many air bubbles on the tank walls. These will reduce the cleaning effect in the beginning phase of ultrasonic cleaning. Turning on the degas function will dissipate the droplets and the air bubbles, usually in 90 seconds, and improve cleaning efficiency.

## 13. On/Off button

After powering on, the LED displays which is the default time (as picture shown), and the unit is ready with the normal working status. Press  button **twice** and the cleaning starts. The cleaning stops when the timer counts down to **00:00**. If the unit needs to be stopped before the timer runs out, press  button **twice**.

## 14. ▼ Decreasing button

- When  is used, each pressing reduces the timer by 1 minute.
- When  is used, each pressing reduces the set temperature down to the next level.
- When  is used, each pressing reduces the solution set time by 20 minutes.


## 15. Solution usage timer set button

Each time that  button is pressed, the timer is increased by 20 minutes. For example, 1:20 indicates that the solution is set to be used for 1 hour and 20 minutes. The maximum time can be set is 9 hour 60 minutes.

## 16. Timer quick set button


Press  button, LED display shows **05:00**. Each pressing increases the time by 5 minutes.

## OPERATION

1

Remove the cover and put it on the cover holder vertically. Put items in the basket, then put them in the cleaning tank. Add water to a level between MIN and MAX and above the area to be cleaned. If the unit is turned on without water, ultrasonic energy will not be absorbed. This model also has high power. Once on for over 15 seconds, it may damage the unit or severely reduce the life of the unit.


2

Connect the power cord (L) and turn the switch on (I). LED display shows **05:00**. This is the most common timer setting and working state. If the timer needs to be adjusted, press  button. The timer can be set between 1 to 30 minutes. Common cleaning times are 5, 10 or 15 minutes.

Prolonged cleaning time may result in:

- Loosening of the screws if used.
- Increasing pre-existing cracks.
- Peeling of coating which was already separated.

3


Press  **twice** to start cleaning.

During cleaning, buzzing sound can be heard from the cleaning tank, indicating that the cleaning is underway. Closing the lid will reduce the noise level.

The digital timer will count down to show the remaining cleaning time. When it displays **00:00**, the cleaning is done.

To stop cleaning at any time, press  button **twice**.

The unit is designed with overheating protection to avoid extended and continuous operation for too long. When the red light is on, the unit will not turn on even if it is forced to. The unit will be idled for about 20 minutes and can be restarted after the green light comes on.


4

When the heater is hended, press  button to select one of the 5 preset temperatures. Press  button **twice**, the icon  flashes, indicating the heater is turned on. When the water reaches the set temperature,  dims, indicating the heater is turned off.

To ensure safety, the unit is designed with dual protections. The heater will be turned off automatically if it has been running for 50 minutes.

If solution is added, to speed up the mixing and to use less solution, use degas function. Press  button to set the timer to 10 minute. Then press **Function** button . The icon  illuminates. Press  **twice** to start degassing. It changes over to normal cleaning automatically after 90 seconds.

5

When debris “smoke” is no longer visible, cleaning is done. If additional cleaning is needed, reset the timer and repeat the steps above.

When cleaning is completed, switch off the power, open the lid and retrieve the basket and the items.

Alternatively, raise and suspend the basket above water to drain water.

Connect the hose to the drainage pipe as shown.


Open the drainage valve to drain the dirty water.

Clean and wipe the tank dry.

Close the drainage valve.

## MAINTENANCE AND CLEANING


In order to ensure an optimum performance and reliability of the unit, we recommend the following:

### **Do not turn on the unit without water in the tank.**


Even though the unit is designed with multiple protections, if it is turned on for over 30 seconds without water in the tank, it may damage the unit or severely reduce the life of the unit. If heating over 3 minutes, it will damage the heater.

This unit is equipped with a thermal cutoff. If the cleaning function is turned on without water in the tank, the temperature will rise continuously. The thermal cutoff will be actuated after approximately 15 minutes respectively to ensure safety. The unit will be disabled until it's repaired.


### **Do not run the unit for extended time or continuously.**

The unit is designed with overheat protection. In the unit has been running for 45 minutes, it is recommended to stop the unit for about 20 minutes to prolong the life of the unit.


### **Do not keep water in the cleaning tank for a long time.**

After Cleaning is completed, open the drainage valve to drain the dirty water. Clean and wipe the tank dry.

### **Do not spray water over the housing and keyboard.**

Use a towel to wipe and dry the tank and the housing.

### **Do not expose the unit under direct sunshine for a long time.**

Keep the unit in a dry, cool and ventilated area.

- Dry the tank and unit with a towel. Do not wash the unit completely with water to prevent the risk of shock.
- Do not load the tank excessively; it could reduce the efficacy of the bath and damage the unit.
- Do not put very heavy or large objects in the bottom of the tank.
- Use chemical cleaning products compatible with the objects to be treated. Avoid corrosive solvents.
- When items (especially valuable items) are being cleaned for the first time, carry out a trial run first with just one item. Ask the manufacturer what cleaning solution is recommended.
- During normal operation the bath and the base of the machine can heat up considerably. Temperatures up to 70/80 °C are completely normal, and can be reached after a long and continuous operation with the lid on. If you wish to reduce the temperature, remove the lid and let the machine rest for a few minutes before starting a new cycle.
- After using the unit, clean the tank and the external body with a damp cloth. Always unplug the unit before cleaning it.
- Keep the machine in a cool and dry place

## PRECAUTIONS

**IMPORTANT- The input voltage of the Ultrasonic Bath transducers can reach hundreds of volts. Therefore, it is very important to observe the following Safety Regulations:**

- 1) Do not open the lid when the power cord is connected to the power supply.
- 2) Do not spill water or any other conducting fluid over the unit's housing.
- 3) Always connect the unit to grounded power outlets.

- ◆ Connect the unit to a grounded 230 V, 50/60 Hz power outlet.
- ◆ Do not use inflammable solvents as cleaning solutions.
- ◆ Keep the unit unplugged before filling it. Do not fill the tank above the *max.* mark. Do not start up the unit when the water level is below the *min.* mark.
- ◆ Do not allow children or unskilled people to handle the unit.

- ◆ Never submerge the unit in any liquid to clean it.
- ◆ Ultrasonic cleaning is not suitable for organic pearls.
- ◆ Do not open the unit while it is connected to the power supply.
- ◆ The machine lid may be hot after a cycle. Proceed with caution when removing it.
- ◆ If you are not sure of the effects of the cleaning solution on a given type of item, carry out a previous test with a single item.
- ◆ Do not attempt to remove any screw or panel. The user does not need to manipulate anything inside.
- ◆ To prevent electric shock, never touch the power cord or socket with wet hands.
- ◆ Never operate the unit if the plug or the power cord is damaged. If the unit does not work properly, is damaged or has fallen into the water, refer it to a technical service for examination and repair.
- ◆ Do not touch the unit if it has fallen into the water. Unplug it immediately.
- ◆ Use the unit properly according to the user's guide. Do not use accessories not recommended by the manufacturer.
- ◆ The unit should never be left unattended when it is plugged.
- ◆ Do not operate the unit next to heat sources such as radiators.
- ◆ Install the unit on a dry and flat surface.



## CLEANING SOLUTIONS

Cleaning solutions can be classified into two types: those with water as the basic element and cleaning solution, and solvents.

The most common cleaning method is with water, as it has an acceptable capacity to dissolve inorganic substances. However, it is extremely limited in dissolving organic substances. For this type of compound, water is a good medium for the dispersion of waste particles, and a good support for specific detergents which dissolve organic substances.

In most cases a mixture of water and domestic washing up liquid can be used as a cleaning solution. Nevertheless, when the quantity or nature of dirt is resistant to this type of cleaning solutions, specific solvents for ultrasonic baths have to be used. **Under no circumstances should inflammable solvents or environmentally unfriendly products be used.**

## PRINCIPLES OF ULTRASONIC CLEANING

Millions of tiny air bubbles are generated within liquid by high frequency vibration.

The air bubbles burst when in contact with object.

The impact dislodges the debris to achieve the cleaning effect.

## TECHNICAL SPECIFICATIONS

Power supply:	230 V, 50 Hz
Power:	300 W
Frequency:	35 kHz
Outside dimensions:	515 x 305 x 240 mm
Tank dimensions:	323 x 205 x 99 mm
Tank capacity:	6 L <i>min: 3.35 L / max: 5 L</i>
Max. item length:	340 mm
Tank material:	Stainless steel
Timer:	1 - 30 minutes
Weight:	7.5 kg

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*TALLERES MESTRAITUA S.L. is a company manufacturing tools and appliances for the dental prosthesis sector, from the business address shown at the bottom of the page.*

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**I HEREBY STATE:** That the **ULTRASONIC CLEANING BATH R-100136**, with the following characteristics:

Tensión de alimentación / Voltage:	230 V, 50 Hz
Potencia / Power:	300 W
Frecuencia / Frequency:	35 kHz
Dimensiones exteriores / Outside dimensions:	515 x 305 x 240 mm
Dimensiones de la cuba / Tank dimensions:	323 x 205 x 99 mm
Capacidad de la cuba / Tank capacity:	6 L min: 3.35 L / max: 5 L
Longitud max. objeto / Max. item length:	340 mm
Material de la cuba / Tank material:	Stainless steel
Temporizador / Timer:	1 - 30 minutes
Peso / Weight:	7.5 kg

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*Meets the low voltage requirements set out in Directive 2014/35/UE and in Directive 2014/30/UE concerning electromagnetic compatibility.*

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*In witness whereof, the company issues this declaration of conformity in Sondika, on 11 December 2017.*

Firma / Signed:

Ignacio Mestraitua,  
Director Gerente / Managing Director

**MESTRA<sup>®</sup>**

TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tel. +34 944530388 - Fax +34 944711725

mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA