ACRYORT SC

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 30°C
88°F
MAX
30°C Limite superiore di temperatura

 Fabbricante

(GB) SELF CURING ACRYLIC RESIN FOR ORTHODONTIC PLATES AND BITES.

DIRECTIONS:

Self curing acrylic resin for orthodontic appliances and bites. It is also particularly suitable for the spray-on technique.

CONTRAINDICATIONS:

- Direct contact of non-polymerised material in the intraoral area.
- In case of proven allergy to components of Acry Ort Sc.
- All uses not indicated.

PREPARATION: It is recommended to work at an ambient temperature of 22°/23°C.

A) For the doughing technique: mix in a ratio of three parts powder to one part liquid (2:1 by weight).

Stir vigorously for about 20 seconds and leave to rest for about 2 minutes. Processing time is around 5-6 minutes, after which the appliance is immersed in a pot or polymeriser at a pressure of 1.5-2 bar for 10 minutes at a temperature of 40°/45°C. If the product is not subjected to polymerisation under pressure, after 9 minutes an exothermic reaction will start until total polymerisation in a free environment.

B) For the direct spray-on technique: submerge the plaster model in water until completely saturated to avoid porosity in the finished product.

Air dry and insulate. After preparing the model to be treated (fastening of hooks and arc screws with hard wax), spray the polymer directly on the model using the container. Pour a sufficient quantity of the monomer to obtain a completely soaked product. Repeat the operation several times until the product is completed. After modelling the device with the resin, it is recommended to remove any excess using a sharp spatula.

Processing time for this stage is around 5 minutes. Immerse the device in a pot or polymeriser at a pressure of 1.5-2 bar for 10 minutes at a temperature of 40°/45°C. If the product is not subjected to polymerisation under pressure, after 9 minutes an exothermic reaction will start until total polymerisation in a free environment.

FINISHING AND POLISHING: Finish as usual, with burs, rubber points, brushes and with pumice.

Mirror polish with the help of brushes and polishing products (Acry Lux Crema, Acry lux Liquido, Acry Lux Pasta).

N.B.: Observing the above procedures will result in a denture with high stability and excellent tensile strength.

Powder: Store in a cool and dry place.

Liquid: Store at a maximum temperature of 30°C.

The product is designed for professional use in the dental sector. Follow the instructions for use. The manufacturer shall accept no responsibility for damages arising from different or inadequate use. The user must personally check the suitability of the product for the intended use, particularly if these uses are not reported in the instructions for use. Do not use after expiration.

Warning: Any serious accident that may occur with the device must be reported by the user to the manufacturer and the competent authority of the Member State where the user and/or patient are established.

LIQUID: Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid breathing vapours. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. Store in a well-ventilated place. Keep cool. Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains methyl methacrylate, N-N-dimethylparatoluidine.

POWDER: Do not breathe dust. IF ON SKIN: Wash with plenty of soap and water.

Contains benzoyl peroxide, methyl methacrylate. May produce an allergic reaction.

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Ruthinium® - Dental Manufacturing S.p.A.