VOCO

Ufi Gel P

Instructions for use

Product description:

Ufi Gel P is an A-silicone based permanently soft relining material for hand mixing. Its autopolymerising, cold-curing properties assure that **Ufi Gel P** is simple and quick to use for direct as well as indirect relining. The set contains an adhesive liquid to achieve a firm bond between denture and relining material, the relining material itself, as well as glazing base and catalyst for smoothing rough surfaces. **Ufi Gel P** is suitable for all PMMA-based dentures. A freshly made **Ufi Gel P** relining can be repaired and corrected.

Fields of application:

- Permanently soft relining for total and partial dentures
- to relieve pressure from pressure sores
- to dam the palatal vibrating line
- in cases of flabby ridge and insufficient adhesion
- to cushion sharp-edged alveolar processes
- to support the healing process in implantology

Method of application:

1. Preparation and coating of the denture:

Clean denture thoroughly with a brush and dry. Always remove any existing previous relinings. Grind down all areas to be relined by at least 1 mm to allow a later relining layer of at least 2 mm thickness. Prepare vestibular, lingual and dorsal rims, e. g., in the form of a chamfer. Round off the rims of the denture. Clean denture and dry by air. Do not use water or a steam jet. Then degrease **exclusively with Isopropanol or pure alcohol (minimum 90 %)** and aerate for 1 min. Humidity, additives (such as those in cavity cleaning agents), saliva, monomer-containing synthetic liquids or touching will diminish the adhesion of **Ufi Gel P** to the denture.

Prior to first use, exchange the bottle cap of the adhesive against the cap with the brush.

Do not use the brush cap for the alcohol.

Apply the adhesive to the surfaces to thus prepared to which **Ufi Gel P** is meant to adhere. The surface has to be covered completely and evenly with the adhesive.

Attention: Close bottle immediately and store in an upright position. Aerate denture for 1 min. approx.

Dispense **Ufi Gel P** base and catalyst in ratio 1 : 1 from the tubes (equal length, equal weight) and mix within 30 s to a homogeneous and bubble-free consistency by spreading it evenly.

Close the tubes immediately and take care not to mix up the tops:

catalyst	=	light blue
base	=	red

Apply mixed material first to the lingual and dorsal, then to the vestibular rim. Then coat the denture base. Ensure that all areas carrying adhesive are completely covered with silicone. On areas to be relined, **Ufi Gel P** should be at least 2 mm in thickness. **Ufi Gel P** has a working time of 2 min after commencement of mixing.

2a. Direct relining in the mouth:

1 min after coating with **Ufi Gel P**, the denture can be inserted into the patient's mouth.

Have the patient gently close the mouth in normal occlusion (1 min). Then the patient should carry out masticatory and swallowing movements (5 min).

After 6 min in total the denture can be removed from the mouth.

Important: Observe the above time intervals strictly.

Finishing should be carried out 10 min after removal of the relined denture from the mouth. Any earlier finishing has to be carried out with extreme care.

2b. Indirect relining in the laboratory:

When relining has to be carried out in the laboratory, isolate the moistened plaster model with a formaldehyde-free alginate-based separating agent for synthetic materials (not recommended: silicone-based separating agent) according to the instructions for use of the separating agent.

Place the denture coated with **Ufi Gel P** on the model and fix both in the relining device or casting flask. Shape functional rims with a metal spatula or a clean plaster knife. After shaping, let **Ufi Gel P** set for 2 - 5 min, then place the relining device containing the denture and the model in the pressure unit for 15 min at 40 °C - 45 °C.

Afterwards, leave the relining device to cool for 20 min.

Attention: Denture, plaster model or relining device must not be worked on with a steam jet during any of these phases.

3. Finishing the relining:

Any excess material can be removed with fine sharp (cuticle) scissors or a scalpel. Finishing of the transition area: If necessary, carry out grinding after 30 min, using polishing disks (REF 2049) and/or abrasive strips. Any earlier grinding process needs to be carried out with great caution. The **Ufi Gel P** set contains 1 polishing disk REF 2049. Abrasive strips are available in varying grain sizes. The most suitable sizes are 80 (coarse) and 120 (medium). Use medium to low rotational speed (5,000 - 8,000 RPM) and apply only very light pressure when finishing the transition area. This prevents heat generation in the relining as well as excessive abrasion of the denture material.

Then clean with alcohol, allow any warmed silicone material to cool down and apply glazing.

4. Corrections:

Immediate correction of the relining is possible. Rinse the relining with water and dry, then apply fresh material, shape it and let it set in the mouth for 7 min or in the pressure unit for 5 min.

Any later correction (within two weeks) of the relining can be carried out after removal of a top layer of **Ufi Gel P**.

5. Sealing:

To seal and smooth the treated surfaces use glazing.

For better dosing, hold the bottle vertically with the dropper down and let air bubbles rise to the top. Homogeneously mix equal numbers of drops of glazing base and catalyst with a brush on a mixing salp. Thinly apply the mixture with the brush on the appropriate margins, i.e. on the dry, freshly cut and ground areas. **Ufi Gel P** glazing remains workable for 1 min 30 s. The setting time at room temperature is 10 min. Cooler conditions prolong the setting time. The setting time can be shortened to 3 min approx. by using a pressure pot (40 °C - 45 °C).

Indications, precautionary measures:

- To prepare the denture and improve the adhesion, additional roughening of the denture material can be done by sandblasting. Afterwards wipe with pure alcohol.
- Storage of **Ufi Gel P** in a refrigerator prolongs the working and setting times.
- Strictly observe the mixing ratio stated. The working and setting times cannot be modified by using smaller or larger amounts of the catalyst paste.
- To achieve optimal adhesion, always comply with the setting times until setting is complete.
- Seal only rough, finished areas (functional rim) with glazing to prevent accumulation of food residues or discolouration.
- Close glazing bottles immediately after use:


glazing base	=	white cap
glazing catalyst	=	red cap

- For finishing, use the enclosed trimming disc, REF 2049, or other suitable burs or polishing discs for flexible materials. Do not inhale any silicone or synthetic dust when finishing the denture or the relining. Wear a protective mask, if necessary.
- **Ufi Gel P** contains modified polydimethylsiloxanes and platinum catalyst, the adhesive contains butanone and methacrylates. **Ufi Gel P** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- Cleaning of the relined denture should be carried out regularly with a special soft denture brush or toothbrush and a non-abrasive cleanser (e.g. neutral soap). Then rinse thoroughly.

Storage:

Store at temperatures of 4 °C to 23 °C. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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