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1. DEFINITION

Ventura flow is a flowable, light-curing, radio-opaque microhybrid composite with extremely low shrinkage stress. The composite material meets the EN ISO 4049 requirements.

2. COMPOSITION

Inorganic filler (barium glass 53 % wt, 28 % by volume, particle size between 50 – 3000 nm), BisGMA, BisEMA, catalysts, inhibitors, pigments.

3. INDICATIONS

1. Restoration of class I, II, III, IV and V cavities on anterior and posterior teeth
2. Extraorally tempered inlays and onlays (direct and indirect), veneers
3. Aesthetic corrections of interdental spaces, enamel hypoplasia, discolorations etc.
4. Blockings, fractures.

4. CONTRAINDICATIONS

Opened pulp, pulpitis, known allergy to methacrylate.

5. SIDE EFFECTS

In individual cases, contact allergies have been described in the case of products of a similar composition. To avoid pulp reactions, it is advisable to cover the exposed dentine in the pulp area with suitable underfilling material (preferably with calcium hydroxide). Note: **ventura flow** does not contain TEGDMA or HEMA.

6. INCOMPATIBILITY WITH OTHER MATERIALS

Avoid underfilling materials which may hinder polymerisation owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol, fall under this category.

7. PROCESSING STAGES

7.1. Tooth cleaning

Brush the tooth to be treated and the teeth neighbouring it with a fluoride-free toothpaste. Clean interdental spaces using strips and dental floss if necessary.

7.2. Shade selection

Under natural light, determine the shade using the VITA(R) SHADE GUIDE.

7.3. Drying

Dry sufficiently, apply rubber dam.

7.4. Cavity preparation

Prepare the cavity in the usual way. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity.

7.5. Underfilling

In the case of preparations near pulp, it is advisable to protect the pulp using an underfilling (e.g. calcium hydroxide). Do not use zinc oxide eugenol!

7.6. Etching

Apply an etching gel (e.g. **ventura etching gel**) to the enamel and dentine according to the manufacturer instructions.

7.7 Adhesion

Apply a light cure adhesive (e.g. **ventura unionbond 2**) according to the manufacturer instructions.

8. RESTORATION WITH VENTURA FLOW

Gradually adapt Ventura Flow up to a maximum layer thickness of 2 mm. Light-cure every layer for 40 seconds. When using matrices, after removing them expose filling to light again from lingual and buccal. **Ventura flow** can be used to care for small cavities or as a liner. Note: **ventura flow** produces an extremely low shrinkage stress and shrinkage.

Warning: After each use, turn the syringe screw back by one turn and re-attach the screw cap. Protect light-curing products from strong sources of light. The specified exposure times refer to halogen or LED light-curing devices with a minimum light intensity of 500 mW/cm² and a wavelength of 400 – 500 nm. They apply to a maximum layer thickness of 2 mm and to all shades. The necessary exposure time may vary depending on the light source and its instructions for use. If in doubt, check the lamp's light output and the necessary exposure time before operation in vitro. During polymerisation, an inhibition layer, which must not be touched or removed if other composite layers are to be applied, forms on the surface.

Light power	≥500 mW/cm ²	≥800 mW/cm ²
Color opaco	40 sec	30 sec
All other colours	40 sec	20 sec

9. FINISHING

Trim filling with 40 µm and 12 µm diamonds, points or burs. Polish to high gloss with discs or silicone points. Interproximal finishing is accomplished by fine grit finishing strips.

10. STORAGE

Protect light-curing products from strong sources of light and heat! In the case of room temperatures above 28° C, it is advisable to store the products in the refrigerator. Prolonged temperatures above 28° C can shorten the shelf life of the product. Under these conditions, the shelf life is 48 months. The expiry date is specified on the packaging.

11. PRECAUTIONARY MEASURES

Close containers after each use using the right lid. Keep out of reach of children. When using tips, ensure before use that the cap has been removed. Commercially available medical gloves do not provide protection against the sensitisation effect of methacrylates. If the product comes into contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice.

12. EMERGENCY PROCEDURES

Contact with the skin: wash with soap and water. Contact with the oral mucous membrane: rinse with large quantities of water. Contact with the eyes: wash immediately with large quantities of water and consult an ophthalmologist.

13. HYGIENE

Use the application instruments for only one patient. Clean and disinfect all other instruments after each use. To avoid contamination, dispense the materials at a distance from the patient.

14. WARRANTY

Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. We accept no liability for other damage, namely that caused by non-compliance with the instructions for use or other improper handling or unintended use of a product. Before using the products, the user should check whether they are suitable for the intended purpose. He alone assumes all of the risks associated with using the product and bears sole responsibility for any damage caused as a result of this.

Instrucciones for use issued in February 2018 (2018-02).

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