ventura Impress 2 Light Body

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EN INSTRUCTIONS FOR USE

Ventura Impress 2 is a system of materials for high precision dental impressions based on addition curing silicones (polyvinylsiloxane), particularly indicated for the double impression technique. Both the moldable putty, **Ventura Impress 2 Putty Soft**, as well as the fluid paste **Ventura Impress 2 Light Body**, packaged in cartridges, possess, once polymerized, excellent mechanical properties and a great capacity for reproducing details. Handling is quite simple, since in all cases the base putty and the catalyst have the same density, and the mixture is always made with a 1:1 ratio (by weight or by volume). The different colors of the two pastes provide good control of the mixture. The use of the fluid paste in cartridges greatly simplifies the procedure, producing homogeneous mixtures with a total absence of bubbles and a notable economy in material usage. The hydrocompatibility of **Ventura Impress 2 Light Body** allows it to enter into close contact with humid surfaces without absorbing liquid, making possible a more perfect reproduction of details, particularly on the edges of the preparation. This property also makes possible a better gypsum flow over the surface when preparing the model, penetrating in even the smallest cracks. The polymerization of the polyvinylsiloxane does not produce volatile secondary components. This gives it excellent dimensional stability, which permits the preparation of the model even several days after taking the impression, without clinically notable alterations.

PREPARATION OF THE DISPENSING GUN

1. Load the cartridge by sliding it vertically through the leads provided for this purpose.
2. Press the lever until the plunger rests on the cartridge.

MIXING

1. Unscrew the cartridge cap by turning it 90° (1/4 of a revolution), remove it and discard it, as it cannot be reused.
2. Extrude impression material until it exits both chambers at the same rate.
3. Insert the mixing tip, blocking it by means of a 90° turn (1/4 of a revolution). If an intra-oral injection tip is also necessary, place it over the mixing tip and then place both tips on the cartridge.
4. Lightly press the lever until the mixed material starts to come out of the mixing tip.
5. Once you have finished the procedure, leave the mixing tip on the cartridge until the next use as it serves as a cap.
6. Between uses, store the dispensing gun with the mixing tip facing down.
7. For the next mix using the same cartridge, remove the used mixing tip by turning it 90° (1/4 of a revolution), insert a new mixing tip and proceed.

MAKING THE IMPRESSION

1. The material should be placed in the mouth within 30 seconds from starting the mix, and kept in position for at least 1 minute and 30 seconds (fast set) or 3 minutes (regular set).
2. To avoid distortions in the corrective or second impression, when the tray is inserted, pressure should be applied for no more than 3 seconds, and removal of the impression from the mouth should be done quickly, avoiding oscillating movements.

MODELS

Models can be poured after 30 minutes, and within seven days of polymerization. The impression should be kept at room temperature until the moment of removal.

Ventura Impress 2 Light Body is compatible with all quality gypsums, such as **Fuji Rock, Whip Mix or Ventura**.

ELECTROPLATING (GALVANIZATION)

Ventura Impress 2 Light Body impressions can be galvanized by acid baths of copper or silver alkalines.

PRECAUTIONS

In sensitive subjects, vinyl-poly-siloxanes may cause irritation or other allergic reactions. Contact with latex gloves may prevent vinyl-poly-siloxane from setting. Avoid direct or indirect contact with latex gloves. Instead, wash your hands thoroughly to eliminate all traces of impurities or use vinyl gloves. Silicone adheres to clothes from where it may be very difficult to remove, even by using chemical agents.

STORAGE

In order for the product to maintain its optimum characteristics, it is advisable to store it at temperatures between 10°C and 23 °C and use it before the expiry date.

| TECHNICAL CHARACTERISTICS | Light Body Regular Set | Light Body Fast Set |
|--|--------------------------|--------------------------|
| Type | : Fluid silicone | : Fluid silicone |
| Classification ISO 4823 | : Low viscosity - Type 3 | : Low viscosity - Type 3 |
| Color | : Green | : Orange |
| Mixing ratio base : catalyst | : 1 : 1 | : 1 : 1 |
| Working time at 23 °C | : 2 m 15 s | : 1 m 30 s |
| Setting time | : 3 m 30 s | : 2 m |
| Compression deformation | : approx.3 % | : approx.3 % |
| Recovery following deformation | : 99.8 % | : 99.8 % |
| Linear dimensional variation at 24 hours | : < 0.20 % | : < 0.20 % |
| Shore A Hardness at 24 hours | : approx. 53 | : approx. 48 |
| Instructions for use issued in June 2018 (2018-06) | | |

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MADSPA, S.A. - Río Jarama, 120 - 45007 TOLEDO - ESPAÑA (SPAIN)