ventura similux

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1. DEFINICIÓN

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EN INSTRUCTIONS FOR USE

1. DEFINITION

Ventura Similux is a universal light-curing microhybrid composite for ante-rior and posterior restorations. It is based on BIS-GMA-resin and inorganic filler particles of 0,05-1,5 µm. The total filler load is 81,00 %, the total filler volume is 65,00 %. The composite material meets the EN ISO 4049 requirements.

2. COMPOSITION

Ventura Similux contains: Resin based on BIS-GMA, radio-opaque glass filler, silica, polymer powder, catalyst, stabilizers and pigments.

3. INDICATIONS

Suitable for all cavities.

4. COSMETICAL AND PHYSICAL PROPERTIES

Very aesthetic.
Very high strength.
Abrasion resistance.
Highly radiopaque.
Excellent colour stability.
High gloss polishability.
Highly filled.
Fluoride content.

5. SIDE EFFECTS

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of hard setting calciumhydroxid material.

6. INCOMPATIBILITY WITH OTHER MATERIALS

Do not use in combination with substances containing eugenol because eugenol can impair the hardening of the composite and causes discoloration. Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol.

7. PREPARATION OF TOOTH AND COLOUR MATCHING

Clean the tooth with a fluoride-free polishing paste prior to preparation and colour matching.
Ascertain the tooth shade while teeth are still moist and select the appropriate **Ventura Similux** according to VITA® -Shades.

8. APPLICATION

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction. Best pulp protection in case of deep excavations is the photoactivated disinfection (PAD). If this is not possible, deep excavations should be covered with a thin layer of calciumhydroxide material. Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with **Ventura Etching Gel** and subsequent application of a light cure bonding agent (e.g. **Ventura UniBond 2**) is recommended. Otherwise use a self etching bonding agent before before application of **Ventura Similux**.

9. PLACEMENT OF VENTURA SIMILUX

Apply **Ventura Similux** in the selected shade. **Ventura Similux** is easy to model. Place it by using a non-metallic instrument. Transparent matrix strips may be used.

10. CURING OF VENTURA SIMILUX

Use a light source that emits visible light with a minimum intensity of 1000 mW/cm², holding the tip of the light guide as close to the restoration as possible.
Cure for 20 seconds in case of lighter shades (e.g. A1, A2 and C2).
Cure for 30 seconds in case of darker shades (e.g. A3,5).

For minimizing the effects of polymerization shrinkage, curing in layers is recommended, if the restoration is more than 2 mm deep.

11. FINISHING

Contour with fine diamonds, points or burs. Polish to high gloss with discs or silicone points. Interproximal finishing is accomplished by fine grit finishing strips.

12. MARKING OF THE PACKAGING

The lot number and expiry date are printed on each package.

13. STORAGE

Do not store above 25 °C (78 °F).

14. EMERGENCY PROCEDURES

Contact with the skin: wash with soap and water.
Contact with the oral mucous membrane: rinse with large quantities of water.
Contact with the eyes: wash immediately with large quantities of water and consult an ophthalmologist.

15. HYGIENE

Use the application instruments for only one patient. Clean and disinfect all other instruments after each use. To avoid contamination, dispense the materials at a distance from the patient.

16. ADDITIONAL NOTES

Do not use any resin to adjust viscosity of composite restorative material. Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.
VITA is a registered trade mark of the VITA-Zhanfabrik, Bad Säckinge, Germany.
The material has been developed for the exclusive use of dental professionals.

Instructions for use issued in 2017-12 (December 2017).

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