making work easy









- DE
- EN
- FR
- IT
- ES
- PT
- TR
- RU
- PL
- CS
- ZH
- JA
- KO
- AR
- UK








Silent TS

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














Made in Germany

21-9191 11102019

Silent TS

Nr. 2921-0050 / 2921-1050

DEUTSCH

Inhalt

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ENGLISH

Content

Introduction.....	1
Symbology.....	1
Information for Operators	2
Operating Instructions	
1. Setup and Commissioning.....	2
1.1 Setup.....	2
1.2 Connection to the Extraction Point.....	2
1.3 Electrical Connection	2
1.4 Connecting electrical equipment	2
1.5 External Exhaust Air Route	3
2. Operation	3
2.1 Operating Elements	3
2.2 Switching the Unit ON / OFF.....	3
2.3 Adjusting and Displaying the Extraction Force	3
2.4 Continuous Operation	3
2.5 Automatic Mode	4
2.6 Adjusting the Automatic Start Feature.....	4
2.7 "Change Filter" - Indicator	4
2.7.1 Adjust the sensitivity of the „Replace filter“ indicator.....	4
2.7.2 Activating / deactivating the „Replace filter“ indicator tone	4
3. Cleaning / Maintenance	5
3.1 Seals	5
3.2 Replacing the Dust Bag	5
3.3 Filters	5
3.3.1 Replacing the Fine Particle Filter	5
3.3.2 Replacing the Exhaust Filter	5
3.3.3 Replacing the Electronics Filter.....	6
3.4 Fuses	6
3.5 Self-diagnosis.....	6
3.6 Replacing the Suction Turbine	6
3.7 Replacing the Exhaust Air Labyrinth	7
4. Spare Parts	7
5. Standard Delivery	7
6. Delivery Versions	7
7. Accessories	7
8. Error List	8
Information for Operators	
A. Application Area	10
A.1 Proper Use	10
A.2 Ambient Conditions (in accordance with DIN EN 61010-1).....	10
B. Hazard and Warning Information	10
C. Authorised Individuals.....	11
D. Preparations Prior to Starting	11
D.1 Connecting the Dust Extractor	11
D.2 Connecting a Dust-Generating Device.....	11
E. Repairs	12

F. Disposal Information	12
F.1 Disposing of Consumables	12
F.2 Disposing of the Unit.....	12
F.3 Disposal instructions for countries in the EU..	12
G. Technical Specifications.....	12
H. Liability Exclusion	12
I. Warranty	12

Introduction

We are pleased with your decision to purchase *Silent TS* laboratory dust extractor.

This device sets a new standard with regard to functionality, performance, and ergonomics.

 **Please read the following operating instructions carefully and observe the information they contain in order to ensure a long and trouble-free service life.**

Symbology

The following symbols are employed in these instructions and on the unit itself:

-  **Danger**
This indicates a direct risk of injury.
-  **Electrical current**
This indicates a hazard due to electrical current.
-  **Attention**
Failure to observe the associated information can result in damage to the unit.
-  **Note**
This provides the operator with useful information to make working with the unit easier.
-  **Only intended for indoor use.**
-  **Before opening the unit, disconnect it from the mains power supply by unplugging the power cord from the wall outlet.**
-  **Burn hazard**
Hot surface or objects.
-  **Observe the operating instructions.**

Other symbols are explained as they occur.

Information for Operators

 Using these operating instructions as a starting point, instruct all operators of the unit with regard to the area of application, the possible hazards during operation, and the proper operation of the dust extractor.

Please have these operating instructions readily available for the operators. Additional information can be found in the Section, "Information for Operators", at the end of these instructions.

Operating Instructions

1. Setup and Commissioning

1.1 Setup

Remove the unit and all accessories from the shipping package.

Inspect the delivery for completeness (refer to the "Standard Delivery" Section).

The unit is fully operational upon delivery.

 **Select a setup location for the *Silent TS* where the exhaust air vent is not blocked.**

Where the unit is to be installed in a cabinet, an opening with the following minimum dimensions must be provided:

- Circular opening: Diameter at least 120 mm
- Rectangular opening: At least 170 x 65 mm.

The *Silent TS* is a free-standing unit intended to be set up on the floor (e.g., under a workbench).

 **Particularly after installation of the roller set (refer to the accessories list), the dust extractor may only be operated if it is standing on the floor.**

1.2 Connection to the Extraction Point

Use the included suction hose to connect the unit to the extraction point.

- Push the suction hose onto the extractor connection fitting on the dust extractor (Fig. 1).
- Connect the suction hose to the desired extraction point (e.g., *Dustex master plus*, extractor clamp, etc.).

With the aid of an extractor switch or a Y-junction it is possible to connect up to two suction points, or if necessary, use a universal adapter for connector tubes (see accessories).

 **Never operate the *Silent TS* without a suction hose.**

 **Avoid steep pitches or hanging points along the hose path. Never kink the hose and make sure it is never stretched or under tension when connected.**

1.3 Electrical Connection

 **Before connecting the unit to the wall outlet, make sure the voltage information on the nameplate corresponds to your local power supply.**

Arrange conducting parts (mains outlets, plugs and couplings) and install extension cord such that the protection class (IP) is retained.

- Unroll the power cord.
- Press the OFF switch (Fig. 2A).
- Plug the power cord into the wall outlet (Fig. 3).

The *Silent TS* is now ready for continuous operation.

1.4 Connecting electrical equipment

Electronic equipment (electrical dust-generating equipment) can be connected to the *Silent TS* at the coupler socket at the back of the machine (N, Fig. 4). In the automatic mode, extraction starts automatically when a connected electric device is switched on. Section D.2 at the end of these instructions offers information concerning electrical systems employing different power plug designs.

 **When connecting electrically operated equipment to the dust extractor, please observe the operating instructions and safety information provided with the equipment.**

 **The unit power outlet is to be used only for connecting electric, dust-generating devices to be operated in conjunction with the dust extractor.**

 **The unit power outlet is live once the *Silent TS* has been connected to the power supply – even when the *Silent TS* is switched off. This allows connected devices to be operated without extraction for brief periods.**

 **Power consumption of all devices connected to the unit power outlet may not exceed the rated value printed next to the unit power outlet.**

1.5 External Exhaust Air Route

An external exhaust air route (see accessories) allows the extracted air to leave the laboratory. The installation details are supplied with the external exhaust air route.

 **When the extraction unit is used in conjunction with an external ventilation system, a significant quantity of air is extracted from the room per hour.**


This can create negative pressure within the room which, when using an air dependent naked flame fed by gas, liquid or solid fuel, can cause poisonous gasses (e.g. carbon monoxide) to be drawn into the working area.

It is therefore essential to ensure that the fresh air supply is sufficient and that the environmental air pressure is maintained, this should then be monitored by an authorized specialist (e.g. a certified Gas Service Engineer).

2. Operation

2.1 Operating Elements

Refer also to Fig. 4, 5

-  (A) OFF switch
-  (B) ON switch
-  (C) Automatic mode LED indicator
-  (D) Continuous / automatic operating mode selector switch
-  (E) Continuous mode LED indicator
-  (F) Extraction force / activation sensitivity adjustment switch
-  (G) Extraction power / activation sensitivity display
-  (H) Extraction force / activation sensitivity adjustment switch
-  (K) "Change Filter" LED indicator
- (L) Mains fuse (fig 4)
- (M) Motor protection switch (Fig. 4)
- (N) Unit plug socket (Fig. 4)

2.2 Switching the Unit ON / OFF

The *Silent TS* is switched on at the ON switch (B) (Fig. 2b), and off at the OFF switch (A) (Fig. 2a). Extraction starts and stops in conjunction with the selected operating mode (continuous operation or automatic operation).

 **The equipment's operating state is saved once the suction unit has been operating for longer than approx. 5 sec.**

2.3 Adjusting and Displaying the Extraction Force

The suction performance of the *Silent TS* can be adjusted in 4 stages by the buttons (F) and (H). This permits the unit's performance to be precisely matched to different types of materials being extracted.

The extraction force setting is shown on the display (G).

 Increase the extraction force.

 Decrease the extraction force.

2.4 Continuous Operation


 Continuous operation LED (E)



In the continuous operation mode, the extractor will begin running as soon as it is switched on, regardless of other electrical equipment which may or may not be connected.

The operating mode selector switch (D) is used to switch between the "continuous" and "automatic" operating modes.

The selected operating mode is indicated by LED (C) or LED (E).

Continuous operating mode selection:

- Switch the *Silent TS* on
 - LED  comes on
 - Continuous operation has now been selected
- or

- LED  comes on
 - The automatic mode has been selected.
 - Press the „D“ key
 - LED  comes on
 - Continuous operation has now been selected.
- The extractor will switch on.

The extractor is started and stopped with the ON (B) and OFF (A) switches.

2.5 Automatic Mode

 Automatic Mode, LED (C).


In automatic operation suction only starts when the electric device connected is switched on responsible is operated.

Selecting the automatic mode:

- Switch the *Silent TS* on (Fig. 2b)

- LED  comes on
- The automatic mode has been selected.

or

- LED  comes on
- Continuous operation has now been selected and the extractor starts.
- Press the „D“ key
- LED  comes on
- The automatic mode has been selected.

The extractor starts as soon as you begin operating the dust-generating device.

After the dust-generating device is switched off, the *Silent TS* continues to run for approx. 5 sec., then stops automatically.

2.6 Adjusting the Automatic Start Feature



In the automatic mode, the *Silent TS* is switched on and off by an automatic start system.

This automatic start feature reacts to the current consumed by an electronic equipment plugged into the unit power outlet (N, Fig. 4).

If the current exceeds the activation threshold set at the *Silent TS*, extraction starts. If the current drops below this activation threshold, extraction stops.

i **The factory set parameter for the start-up threshold enables the *Silent TS* to react correctly for most electrical dental equipment. Should this not be the case for a particular piece of equipment (eg. a handpiece), then the start-up threshold can be adjusted accordingly.**

To change the activation threshold:

1. Switch *Silent TS* on.
2. Press and hold  for 3 sec.;
 - LED  and  flash.
 - A “0” flashes on the display.
3. Switch electronic equipment off. For units equipped with a standby mode, switch them to this mode (e.g., for units equipped with handpieces, only switch the controller on without activating the handpiece).
4. Press the  key.
 - A “1” flashes on the display.
5. Switch the electrical appliance on, i.e. activate the hand piece at the speed at which the *Silent TS* is required to start up.
6. Press the  key.
 - An audible signal confirms that the adjustment has been successfully performed.

The *Silent TS* returns to the operating mode (continuous or automatic operation) which was set for the automatic start function prior to the adjustment.

2.7 "Change Filter" - Indicator

LED (K) coming on indicates that either the dust bag or the fine particle filter needs to be replaced.

An acoustic signal will also sound three times in this case. In this case, follow the instructions given in the “Cleaning / Maintenance” section to replace the dust bag or the fine particle filter.

i **If the dust extractor continues to be used after the “Change Filter” indicator goes on, it will automatically stop once the airflow drops below a predefined minimum value. The LED display blinks (K). The indicator can be turned off by switching the *Silent TS* OFF and ON, thus restarting extraction.**



 **Operating the unit with a blocked filtration system can result in hazards and damage to the extractor. Filters must be replaced if the extractor stops due to a full filter.**

2.7.1 ADJUST THE SENSITIVITY OF THE „REPLACE FILTER“ INDICATOR.

i **It may be necessary to reduce responsiveness if the suction support of the connected appliance has a smaller diameter than the *Silent TS*. (<< 32.5 mm).**

The „Replace filter“ sensitivity can be adjusted within limits.


To do this:

- Switch the *Silent TS* off;
- Switch the *Silent TS* on, and in doing so keep the  button firmly pressed until the display begins to blink. (5=factory set parameter)
- Press button ,
 - Reduce level of sensitivity,
 - „Change filter“ signal will appear later.
 - ==> Dust bag will get very full
- Press button ,
 - Increase level of sensitivity,
 - „Change filter“ signal appears more quickly.
 - ==> Dust bag will get less full
- Press button ,
 - Acoustic signal confirms successful adjustment.

2.7.2 ACTIVATING / DEACTIVATING THE „REPLACE FILTER“ INDICATOR TONE

The „Replace filter“ indicator tone can be activated or deactivated.

To do this:

- Switch the *Silent* off;
- While holding down the  and  keys, switch the *Silent* back on and wait for the indicator tone
 - Brief indicator tone: acoustic signal deactivated.
 - Longer indicator tone: acoustic signal active.

3. Cleaning / Maintenance

 Always unplug the unit from the wall outlet before beginning any cleaning or maintenance tasks.

 Examine the power cord regularly, but at least annually, for damage or signs of aging. Immediately replace damaged power cords.

3.1 Seals

To ensure the proper function of the dust extractor, it is vital that the following three seals:

- Profile seal on the dust drawer (Fig. 12a);
- V-ring seal on the dust drawer (Fig. 12b);
- Seal on the upper front panel (Fig. 20);

are not damaged.

Inspect these seals when replacing the associated filters and replace them if they are found to be damaged (refer to the spare parts list).

3.2 Replacing the Dust Bag

The dust extractor may only be operated with a complete filtration system. The dust bag must be immediately replaced with a new one as soon as the “Change Filter” indicator goes on.

 **Failure to replace the dust bag may result in it's rupturing. Depending on the material in the bag, this may lead to hazards for the operator. In addition, the dust extractor may be damaged.**

- Pull the front panel forward and off (Fig. 8).
- Pull the dust drawer out towards the front (Fig. 9).
- Leave the dust bag in the dust drawer and take it to the disposal location.
- Remove the dust bag from the drawer.
- Take off protection film, cover filter opening (Fig. 10).
- Dispose the dust bag accordingly.

 **Always comply with all local ordinances governing proper disposal and accident prevention!**
Depending on the filter contents, protective clothing may be required.

- Insert the new dust bag into the dust drawer, being careful to ensure that the dust bag is properly aligned in the guide grooves and the glue lap faces up (Fig. 11).

 **Use only original Renfert dust bags (refer to the spare parts list).**

- Inspect the dust drawer seal for damage and replace it, as required (Fig. 12a, b).
- Insert the dust drawer, making sure it is straight, and push it fully against its stops. Make sure the dust drawer is seated correctly on the guide (Fig. 13).
- Install the front panel on the bottom hooks first (Fig. 14), then push it into the upper lock bolts, pressing firmly until it locks in place (Fig. 15).

3.3 Filters

 **NEVER operate the dust extractor without the complete filtration system.**

3.3.1 REPLACING THE FINE PARTICLE FILTER

Examine the fine particle filter regularly, but at least annually, and replace it. It must always be replaced if the “Change Filter” indicator (K) continues to go on, despite the dust bag having been replaced.

- Pull the front panel forward and off (Fig. 16).
- Turn the fine particle filter anticlockwise to release it (Fig. 17).
- Pull the fine particle filter straight forwards and out of the unit (Fig. 18) and dispose of it properly.

 **Always comply with all local ordinances governing proper disposal and accident prevention!**
Depending on the filter contents, protective clothing may be required.

- Install a new fine filter and fully insert it over the holder bar (Fig. 19).
- ! **Use only original Renfert fine particle filters (refer to the spare parts list).**
- Turn the fine particle filter clockwise and hand-tighten it in place (Fig. 17).
- Inspect the seal on the front panel for damage and replace it, as required (Fig. 20).
- Snap the front panel in place (Fig. 21, 22).

3.3.2 REPLACING THE EXHAUST FILTER

The exhaust filter primarily traps particles rubbed off the suction turbine graphite brushes. The filter should be replaced annually:

- Press on the four lock tabs on the filter cassette and pull the cassette down and off the unit (Fig. 23).
- Properly dispose of the filter and the filter cassette.

 **Always comply with all local ordinances governing proper disposal and accident prevention!**
Depending on the filter contents, protective clothing may be required.

- Insert a new filter mat in the filter cassette so that the smooth, compressed side of the filter mat faces down or outwards when the filter is installed.
- Install the new filter cassette on the exhaust air vent.
- Make sure the filter cassette is properly seated and locked in place on all four tabs.

 **Use only original Renfert exhaust filters (refer to the spare parts list).**

3.3.3 REPLACING THE ELECTRONICS FILTER

Cooling air for the unit's electronics passes through the electronics filter. The type of material trapped by the filter and the degree of contamination depend on the ambient laboratory conditions.

The filter should be examined annually and replaced, as necessary:

- Press on the two lock tabs on the filter cassette and pull the cassette off towards the rear of the unit (Fig. 24).
- Properly dispose of the filter and the filter cassette.

 **Always comply with all local ordinances governing proper disposal and accident prevention!**
Depending on the filter contents, protective clothing may be required.

- Install a new filter cassette and filter on the opening.
- Make sure the filter cassette is properly seated and locked in place on both tabs.

 **Use only original Renfert electronics filters (refer to the spare parts list).**

3.4 Fuses

Instead, two device overload switches, (L, Fig. 4) are provided to protect these units.

Press the button to reset a switch which has tripped.

3.5 Self-diagnosis

The dust extractor is equipped with a self-diagnosis, which checks various functions and displays errors on the display (G).

 **Before the diagnosis takes place, insert a new dust bag and ensure that the fine filter is clean and the suction pipe is not blocked.**

Activate the self-diagnosis:

- Switch *Silent TS* off;
- Press and hold the operating mode selector switch (D) and switch *Silent TS* on (Fig. 7a).
- For approx. 3 sec.:
 - all LED displays are illuminated,
 - An „8“ is indicated in the display
 - The signalling device gives an acoustic signal.
- A “d” appears on the display during the unit's self-test (Fig. 7b). During this time, the unit's electronics check various internal elements and functions. The suction turbines will briefly run at Level 1 during this test phase.
- If no errors are detected, the *Silent TS* automatically goes into the previously selected operating mode at the end of the self-test (after approx. 10 sec.).
- If an error is detected during the self-test, a flashing “E” will appear on the display, together with a number (1 to 3).

These indicate:

E1: Airflow sensor fault; have the unit repaired.

E2: Plug on grey motor wire out of its socket (Fig. 26) or suction turbine fault. Replace suction turbine if necessary (refer to the spare parts list).

E3: Error in the electronics; have the unit repaired.

In this case, the *Silent TS* remains in the self-diagnosis mode until it is switched off at the OFF switch (A).

3.6 Replacing the Suction Turbine

The suction turbine is encased in an encapsulating housing to form a single unit, thus allowing it to be easily replaced without the need for any tools.

 **Before opening the unit, disconnect it from the mains power supply by unplugging the power cord from the wall outlet.**

 **The motor may be hot. Allow the motor to cool off before proceeding.**

 **The suction turbine may only be operated when it is installed in the unit. Only an authorized electrician may perform a functional inspection of or repair the unit.**

 **Please observe the documents included with the spare part!**

- Turn the lock knob 90° anticlockwise (Fig. 25).
- Remove the motor compartment cover.
- Open the fixing clip and remove the motor plug (Fig. 26).
- Loosen the two internal, grey knurled nuts (Fig. 27) and unscrew them approx. 1 cm. These knurled nuts do not need to be completely removed.
- Release the suction turbine by turning it anticlockwise, then pull it straight out towards the back (Fig. 28).
- Install the new suction turbine (Fig. 28) by pushing it in straight and turning it clockwise to lock it in place (Fig. 29). The connecting wire should be on top (Fig. 28a).
- Tighten the internal knurled nuts (Fig. 27).
- Insert the motor plug and close the fixing clip (Fig. 26).
- Check that the lock knob is in the open position (vertical).
- Place the motor compartment cover on the motor compartment. Two guide pins have been provided on the cover to assist in positioning it correctly.
- Turn the lock knob 90° clockwise (Fig. 31).

 **Always comply with all local ordinances governing proper disposal and accident prevention!**

3.7 Replacing the Exhaust Air Labyrinth

Dust from the suction turbine graphite brushes also settles on the walls of the exhaust air labyrinth. Over time, the labyrinth will become discoloured.

Before opening the unit, disconnect it from the mains power supply by unplugging the power cord from the wall outlet.

Replacing the exhaust air labyrinth:

- Turn the lock knob 90° anticlockwise (Fig. 25).
- Remove the motor compartment cover.
- Pull the exhaust air labyrinth out towards the back.
- Install the new exhaust air labyrinth by pushing it in straight, making sure it is properly seated. The openings in the exhaust air labyrinth must face towards the rear (Fig. 30) and the metal grating downwards.
- Check that the lock knob is in the open position (vertical).
- Place the motor compartment cover on the motor compartment. Two guide pins have been provided on the cover to assist in positioning it correctly.
- Turn the lock knob 90° clockwise (Fig. 31).

Always comply with all local ordinances governing proper disposal and accident prevention!

4. Spare Parts

You can find components subject to wear and the spare parts on the spare part list in the internet at www.renfert.com/p918.

The components excluded from the warranty (such as consumables or parts subject to wear and tear) are marked on the spare part list.

Serial number and date of manufacturing are shown on the type plate of the unit.

5. Standard Delivery

- 1 *Silent TS* laboratory dust extractor
- 1 Operating instructions
- 1 Suction hose
- 1 Dust bag (properly installed in the dust drawer)
- 1 Fine particle filter (properly installed)
- 1 Hose fitting adapter
- 1 Grounded wall plug (No. 2921-0050 only)

6. Delivery Versions

- 2921-0050 *Silent TS*, 230 V, 50/60 Hz
2921-1050 *Silent TS*, 100 - 120 V, 50/60 Hz

7. Accessories

- 2921-0002 Dust bag set (5 bags)
2921-0001 Roller set
2925-0000 Extractor clamp (w./o. glass pane)
2925-1000 Glass pane, w. mounting (for the extractor clamp)
2926-0000 Suction switch, incl. 4 mufflers (w./o. hose)
2921-0003 Muffler set, 2 mufflers
2921-0004 External exhaust air control
2934-0007 90° angled connector *Silent*
15-0823 Suction hose (by the metre)
90003-4305 Hose fitting adapter
90003-4240 Suction hose (3 m), cpl.
90003 4826 Suction hose, antistatic, 3 m, incl. 2 end mufflers
90003-4314 Y-junction
90003-4410 H+Hepa filter *Silent TS*

8. Error List

Error	Possible cause	Corrective action
Dust extractor stops suddenly, for no apparent reason.	<ul style="list-style-type: none"> • Dust extractor continued to be used, despite “Change Filter” indicator. • Unit fuse blown. • Motor protection switch (M, Fig. 4) has been activated due to overload (temperature, electricity). • Motor over-temperature protector switch activated. • Plug on grey motor wire out of its socket or loose. • Suction turbine fault. • In-house fuse blown. 	<ul style="list-style-type: none"> • Check and replace the dust bag and fine particle filter (Chapter 3.2 and 3.3). • Check and reinstall the fuse. • Check the motor protection switch, if necessary reset. • Change the motor if repeatedly activated. • Check dust bag and filter, replace as required. (Sec. 3.2 and 3.3). • Check exhaust air filter, replace as required. (Sec. 3.3.2). • Check whether suction hose is blocked. Remove obstruction as required. • Dust extractor can be restarted after being allowed to cool off for approx. 1 hour. • Check plug (Fig. 26). • Carry out self-diagnosis (chapt. 3.5), at fault „E2“ exchange suction turbine (see spare parts list). • Replace the in-house fuse. • Check the total power consumption of the dust extractor plus dust-generating devices.
Dust extractor stops suddenly, for no apparent reason. E4 is indicated on the display.	<ul style="list-style-type: none"> • The electronics may be overheated due to excessively high ambient temperature, e.g., if the unit is installed in a cabinet. 	<ul style="list-style-type: none"> • Allow the extractor to cool down. • Attempt to improve the ventilation conditions. Refer also to Sec. 1.1.
After switching on the unit, the display immediately indicates E4.	<ul style="list-style-type: none"> • Temperature sensor fault. • After overheating, the device has not yet cooled down sufficiently. 	<ul style="list-style-type: none"> • Have the device repaired. • Disconnect the device and let it continue to cool down.
The self-test returns a suction turbine error, even though the turbine is running.	<ul style="list-style-type: none"> • Filter full. • Suction hose blocked. • Suction turbine fault. 	<ul style="list-style-type: none"> • Check dust bag and fine particle filter, replace as required. (Sec. 3.2 and 3.3.1). • Remove obstruction from the suction hose. • Replace the suction turbine (Sec.3.6).
Insufficient extraction force.	<ul style="list-style-type: none"> • Extraction force level too low. • Filter blocked. • Suction hose blocked. • Dust drawer seal not tight. • Upper front panel fails to close properly. • Suction hose damaged. 	<ul style="list-style-type: none"> • Increase the extraction level (Sec. 2.3). • Check dust bag and fine particle filter, replace as required.(Sec. 3.2 and 3.3.1). • Remove obstruction from the suction hose. • Check the proper seat of the dust drawer (Sec. 3.2). • Check the dust drawer seal for damage and replace as required (Sec. 3.1). • Check the proper seat of the upper front panel (Sec. 3.3.1). • Check the upper front panel seal for damage and replace as required (Sec. 3.1). • Inspect the suction hose and replace as required.

Error	Possible cause	Corrective action
“Change Filter” indicator goes on shortly after filters were replaced.	<ul style="list-style-type: none"> • Fine particle filter blocked. • Suction hose blocked. • Cross-section of the dust-generating device’s hose fitting is too small. • Dust is too fine. • Flow sensor defect. 	<ul style="list-style-type: none"> • Replace the fine particle filter (Sec. 3.3.1). • Remove obstruction from the suction hose. • Adjust the sensitivity of the „Replace filter“ indicator(Chapter 2.7.1.). • Adjust the sensitivity of the „Replace filter“ indicator (Chapter 2.7.1.) • Carry out self-diagnosis (chapt. 3.5), at fault „E1“ equipment requires professional repair.
Dust bag ruptures.	<ul style="list-style-type: none"> • The “Change Filter” indicator’s sensitivity has not been adjusted for the material being extracted. 	<ul style="list-style-type: none"> • Reduce the sensitivity of the „Replace filter“ indicator (Chapter 2.7.1.).
Dust extractor fails to start when a dustgenerating device is being used.	<ul style="list-style-type: none"> • Dust extractor not switched on. • Dust-generating device not plugged into the unit power outlet on the <i>Silent TS</i>. • Wrong operating mode selected. • Automatic start feature’s threshold value too high. • Device fuse has blown. • Automatic start feature fails to detect the attached device’s operation. 	<ul style="list-style-type: none"> • Switch the dust extractor on (Sec. 2.2). • Connect the dust-generating device to the <i>Silent TS</i> (Sec. 1.4). • Select the automatic mode (Sec. 2.5). • Adjust the automatic start feature (Sec. 2.6). • Check and reinstall the fuse.(Sec. 3.4). • Some older handheld instruments do not generate a sufficiently strong signal for the automatic start feature. Switch the dust extractor to the continuous mode.
Dust extractor starts, although a connected device is not being used.	<ul style="list-style-type: none"> • Automatic start feature’s threshold value too low. 	<ul style="list-style-type: none"> • Adjust the automatic start feature (Sec. 2.6).
Dust extractor fails to stop when a connected device is switched off.	<ul style="list-style-type: none"> • Dust-generating device not plugged into the unit power outlet on the <i>Silent TS</i>. • Wrong operating mode selected. • Activation threshold too low. 	<ul style="list-style-type: none"> • Connect the dust-generating device to the <i>Silent TS</i> (Sec. 1.4). • Select the automatic mode (Sec. 2.5). • Adjust the automatic start feature (Sec. 2.6).
No satisfactory setting of the automatic start feature can be adjusted.	<ul style="list-style-type: none"> • Some older handheld instruments do not generate a sufficiently strong signal for the automatic start feature. 	<ul style="list-style-type: none"> • Switch the dust extractor to the continuous mode.
Extraction briefly starts up when the unit is switched on.	<ul style="list-style-type: none"> • Activation sensitivity too low. 	<ul style="list-style-type: none"> • Adjust the automatic start feature (Sec. 2.6).
Panels covering the fine particle filter and the dust drawer are difficult to open.	<ul style="list-style-type: none"> • Lock bolts dirty. 	<ul style="list-style-type: none"> • Clean the lock bolts and lubricate them slightly, if required.

Information for Operators

The following information is intended to assist you, the operator, in safely working with the *Silent TS* in your laboratory.

 **Using these operating instructions as a starting point, instruct all operators of the unit with regard to the area of application, the possible hazards during operation, and the proper operation of the dust extractor.**

Please have these operating instructions readily available for the operators.

A. Application Area

The *Silent TS* is a workbench extractor used to extract dusts such as occur in laboratories, e.g., dental labs.

It is solely intended for commercial use in laboratories and not for private, household use.

The *Silent TS* can be operated both manually as well as in conjunction with connected, electric, dust-generating equipment.

It is possible to connect one or two suction points to the extractor. The second suction point can be connected with a Y-adaptor which is available as an accessory (see chapter „accessories“).

 ***In order to protect the health of persons working in a dental laboratory, the law requires that specialised dust extractors must be used.***

National regulations state the MAK value - maximum concentration of suspended particles in the air.

Find out about the legal national threshold value and the type of dust generated in your laboratory.

A.1 Proper Use

Proper use implies the extraction of dry, non-explosive dusts. Fire-promoting, easily flammable, flammable, or explosive materials may not be extracted with the *Silent TS*.

The extraction of liquids, smouldering, or burning materials is prohibited.

Any other use is not „as intended“. The manufacturer shall not be liable for any damages resulting from any such other use.

The „use as intended“ includes compliance with the operating, servicing and preventive maintenance conditions defined by the manufacturer.

A.2 Ambient Conditions

(in accordance with DIN EN 61010-1)

The unit may only be operated:

- Indoors;
- Up to an altitude of 2,000 m above sea level;
- At an ambient temperature range between 5 - 40 °C [41 - 104 °F] *);

- At a maximum relative humidity of 80 % at 31 °C [87.8 °F], dropping in a linear manner to 50 % relative humidity at 40 °C [104 °F] *);
- With mains power where the voltage fluctuations do not exceed 10 % of the nominal value;
- Under contamination level 2 conditions;
- Under over-voltage category II conditions;

*) Between 5 – 30 °C [41 – 86 °F], the unit can be operated at a relative humidity of up to 80 %. At temperatures between 31 – 40 °C [87.8 – 104 °F], the humidity must decrease proportionally in order to ensure operational readiness (e.g., at 35 °C [95 °F] = 65 % humidity; at 40 °C [104 °F] = 50 % humidity). The unit may not be operated at temperatures above 40 °C [104 °F].

B. Hazard and Warning Information

Only intended for indoor use. The unit is only designed for dry applications and may not be operated or stored outdoors or under wet conditions.

The *Silent TS* laboratory dust extractor is an electric device and, as such, carries with it a certain inherent potential hazard. The unit may not be taken into service until any required alterations to comply with regionally specific power plug configurations have been made. Such alterations may only be performed by a qualified electrician.

The unit may only be operated if the information on the nameplate conforms with the specifications of your local mains power supply.

The mains socket on the unit is only designed for the purposes specified in the Operating Instructions. Connecting other devices may cause material damage. Before connecting another device, turn off both the dust extractor and the other device. Read the operating instructions of the other device and comply with the safety instructions contained in the document.

Please observe the national regulations and permitted exposure to dust in a working environment. Please note EN 60335-2-69, Appendix AA, or initiate appropriate enquiries with your trade association or the responsible authorities.

When extracting hazardous materials, always refer to the relevant safety data sheets.

Always wear protective gear when extracting hazardous materials.

It is necessary to wear suitable, personal protective equipment when emptying the dust bag or cleaning, depending on the type of extracted material.


 Regularly inspect connecting lines and hoses (e.g., the power cord) for damage (e.g., kinks, cracks, porosity) or signs of aging. Units exhibiting damaged connecting lines, hoses, or other defects must be taken out of service immediately.

 Always unplug the unit from the wall outlet before beginning any work on the unit's electrical components.

 Never operate the unit without the complete filtration system (dust bag, fine particle filter, exhaust air filter, electronics filter).

 Make sure the dust drawer is fully closed during extraction.

 Never operate the unit without the suction hose.

Never extract flammable or explosive gases, vapours, dusts.

 Never extract liquids.

 Never extract hot materials.

 The unit may not be operated without supervision.

 All required tools must be appropriately cleaned prior to use when the unit is employed for medical purposes or in conjunction with medicinal materials.

 If the dust extractor is employed to extract hazardous materials, appropriate personal protective gear must be worn and steps must be taken to ensure that the exhaust air is properly ventilated. Please refer to the associated safety data sheets for specific requirements.

 Environmentally hazardous extraction substances. Extracted substances may be a risk to the environment. Dispose of extracted material according to local statutory regulations.

 It is the responsibility of the operator that national regulations during operation and regarding a repeated safety inspection of electrical equipment are complied with. For Germany these are the regulation 3 by DGUV (German Statutory Accident Insurance) in relation with VDE 0701-0702 (Association for Electrical, Electronic and Information Technology).

 Only spare parts and accessories supplied or authorized by Renfert GmbH may be used with this product. If other spare parts or accessories are used, this could have a detrimental effect on the safety of the device, increase the risk of serious injury and lead to damage to the environment or the device itself.

 Information on REACH and SVHC is available on our website www.renfer.com, in the Support area.

C. Authorised Individuals

Only properly trained individuals may operate and service the *Silent TS*.

Minors or pregnant women may only operate and service the *Silent TS* if they are wearing appropriate protective gear, particularly if the unit is being used to extract hazardous materials.

Any repairs not specifically described in these operating instructions, in particular the inspection of the suction turbine outside the unit, may only be performed by a qualified electrician.

D. Preparations Prior to Starting

D.1 Connecting the Dust Extractor

The unit is supplied with a power cord and two-contact plug with a ground contact (either DIN 49441 or NEMA, depending on the particular model in question).

The unit may not be taken into service until any required alterations to comply with regionally specific power plug configurations have been made.

Such alterations may only be performed by a qualified electrician.

The unit may only be connected to wall outlets connected to a ground circuit system.

Before operating the unit, make sure the voltage information on the nameplate corresponds to your local power supply.

D.2 Connecting a Dust-Generating Device

The unit is equipped with a power outlet on the back, designed for a two-contact plug with a ground contact (either in accordance with DIN 49441 or NEMA), to which dust-generating equipment can be connected (N, Fig. 4).

A suitable plug is supplied for units with a nominal voltage of 230V and plug socket according to DIN 49441 (Fig. 32). This connector can be used to make an adapter for your local plug configuration.

 **This adapter may only be made by a qualified electrician!**
The adapter must not interrupt the ground circuit system!

E. Repairs

Repairs may only be performed by qualified electricians or authorized dealers.

Depending on the material trapped by the filters, protective gear may need to be worn during repairs.

F. Disposal Information

F.1 Disposing of Consumables

Full dust bags and filters – including motor and electronics filters – must be disposed of under compliance with locally applicable regulations.

Depending on the material trapped by the filters, protective gear may need to be worn during disposal.

F.2 Disposing of the Unit

The unit must be disposed of by an authorized recycling operation. The selected firm must be informed of all possibly health-hazardous residues in the unit.

F.3 Disposal instructions for countries in the EU

To conserve and protect the environment, prevent environmental pollution and improve the recycling of raw materials, the European Commission adopted a directive that requires the manufacturer to accept the return of electrical and electronic units for proper disposal or recycling.

Within the European Union units with this symbol should not therefore be disposed of in unsorted domestic waste:

For more information regarding proper disposal please apply at your local authorities.

G. Technical Specifications

Mains voltage:	220 - 240 V, 50/60 Hz 100 - 120 V, 50/60 Hz
Dust extractor power consumption:	1400 W (230 V) 1400 W (120 V) 1000 W (100 V)
Unit power outlet maximum connecting value:	2000 W (230 V) 360 W (120 V) 500 W (100 V)
Total connected power:	3400 W (230 V) 1800 W (120 V) 1500 W (100 V)
Mains input fuse (Fig. 4, L):	2 x 15 A (T)
LpA *) (at max. volume flow):	63 dB(A)

Volume flow, max.:	3300 l/min (120 V / 230 V) 2900 l/min (100 V)
Vacuum pressure, max.:	200 hPa [2.9 psi]
Filter surface area, fine particle filter, approx.:	0.8 m ² / [1240 sq inch]
Fill volume, dust bag, approx.:	7.5 l [2 US gal]
Dimensions (width x height x depth):	224 x 591 x 564 mm [8.8 x 23.3 x 22.2 inch]
Weight (empty), approx.:	26 kg [57 lbs]
Ø Suction fittings:	
Internal:	32 mm [1.26 inch]
External:	40 mm [1.57 inch]

*) Sound pressure level in accord. with EN ISO 11202

H. Liability Exclusion

Renfert GmbH shall be absolved from all claims for damages or warranty if:

- **The product is employed for any purposes other than those cited in the operating instructions;**
- **The product is altered in any way other than those alterations described in the operating instructions;**
- **The product is repaired by other than an authorized facility or if any but Renfert OEM parts are employed;**
- **The product continues to be employed, despite obvious safety faults or damage;**
- **The product is subjected to mechanical impacts or is dropped.**

I. Warranty

Provided the unit is properly used, **Renfert warrants all components for 3 years.**

i **The suction motor has a guarantee of 3 years, with a maximum of duration of 1000 operating hours (motor running time).**

Warranty claims may only be made upon presentation of the original sales receipt from the authorized dealer.

Parts which are subject to natural wear and tear (wear parts) and consumables are excluded from the guarantee. These parts are marked in the spare part list.

The warranty is voided in case of improper use; failure to observe the operating, cleaning, maintenance, and connection instructions; in case of independent repairs or repairs by unauthorized personnel; if spare parts from other manufacturers are employed, or; in case of unusual influences or influences not in compliance with the utilization instructions.

Warranty service shall not extend the original warranty.

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Renfert GmbH • Untere Gießwiesen 2 • 78247 Hilzingen/Germany