Ref. 080543

- ① Panel superior / **Top plate**
- ② Cuerpo principal / **Main unit**
- ③ Manómetro / **Pressure gauge**
- ④ Resistencia / **Heater**
- ⑤ Anillo superior / **Upper sliding**
- ⑥ Anillo inferior / **Lower sliding**
- ⑦ Placa inferior / **Bottom plate**
- ⑧ Apoyo modelo / **Bead case**
- ⑨ Columna / **Support fixture**
- ⑩ Carcasa resistencia / **Heater case**
- ⑪ Conexión aire / **Air connection**
- ⑫ Regulador / **Regulator**
- ⑬ Etiqueta / **Label**
- ⑭ Toma de corriente / **Power supply**
- ⑮ Muelle / **Spring**

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The pressure former ref. 080543 is a practical machine with many applications in dental laboratory. For optimum performance and a long-lasting life, please read the following instructions.

Strong pressurization allows precise molding, regardless of model trimming, as well as various thickness (0.5 to 4.0 mm), and size (127, 125, 120 mm) square and round sheet can be used.

## Installation

- Carefully unpack the unit.
- Choose a rigid, flat surface, away from sources of vibration or heat, for the location of the unit. Place the unit at a height that allows you to comfortably handle objects inside the chamber.
- Connect the electrical plug to a 220-240 V, 50/60 Hz power supply with an earth connection.
- Connect air tube (Ø 8) to the air connection (11) on the bottom of the back of the body.
- Connect power cable to the power supply (14) on the bottom of the back of the body.
- Remove the support fixture (9) and fixed cable tie.
- Rotate the lowering sliding (6) to the lift (90°) to look to the front.

## Operation

- **Heater:** Sheet heating and softening.
- **Cylinder:** Up & down operation and vacuum.
- **Pressure:** Pressure ejection and model moulding.

## A. Preparation before use

- ① Connect the power cable to the main unit and the outlet.
- ② Connect air tube (Ø8) to air connection of the main unit. (Check pressure gauge operation)
- ③ The factory setting is 5 bar (If the pressure gauge does not work, check the air compressor in use.)
- ④ Rotate the lower sliding in the front direction and place the bead case on the bottom plate of the main unit.

## B. How to use

- ① Select the sheet that best suits your purpose. (Temporary, Splint, retainer, bleaching, and so on)
- ② Place the model you want to make in the middle of the bead case.
- ③ Place a round ○ or square □ type sheet on the lower sliding to fit the guideline, and then put the upper sliding on.
- ④ Press the heater switch on the front top plate. (The lamp is turned on and the heater is in standby mode.)
- ⑤ Rotate the lower sliding to the left (90°). At this time, the heater automatically turns on.
- ⑥ Wait for the sheet to hit the model height or halfway. (Time range is between 50 to 90 seconds. The operating time differs depending on the material and thickness.)
- ⑦ Rotate the lower sliding to the right (180°). (Model location in the center of the body)
- ⑧ Press the cylinder switch on the front top plate. (The lamp is turned on and the air case comes down and presses the upper and lower sliding to create a vacuum standby state.)
- ⑨ Press the pressure switch on the front top plate. At this time, the air is injected into the air case and the sheet is completely molded into the model. (The operating time is about 3 to 5 minutes depending on the material and thickness.)

- ⑩ Turn off the front pressure switch and turn off the cylinder switch on the front top plate. (The sheet is fully seated on the model and comes up.)
- ⑪ Rotate the lower sliding to the left (90°).
- ⑫ Remove the upper sliding and take out the molded model.
- ⑬ Turn off the heater switch on the front plate.

## C. Management and storage methods

- ① Disconnect the power cable from the outlet.
- ② After use, check the equipment is clean.
- ③ Do not place anything on the heater.

## Maintenance and storage

- ① Please keep out of reach of children.
- ② Please use in a place separated from fire and water.
- ③ Temperature: -10~70 °C
- ④ Humidity: 10~100 %
- ⑤ Pressure: 50~106 kPa

## Precautionary measures

- ① There is a risk of burns because the heater area is heated after using this product. Do not touch heater area after operation.
- ② Do not disassemble or repair this product. To avoid electric shock, do not open the case. Refer servicing to qualified personnel only.
- ③ After using this product, must remove the power cable from outlet
- ④ Contact with sheet softening may cause sheet deformation and burns.
- ⑤ Do not use for other purposes.
- ⑥ Avoid strong rays, although this vacuum former has been designed pursuant to international EMI standards.


## Cleaning

Clean the product and each part with a soft cloth and air.

## Components


## Symbols

## Fault coping method

- ① **Power failure:** Check fuse and power switch.
- ② **Heater failure:** Check heater switch, limit switch, heater.
- ③ **Pressure failure:** Check compressor, pressure switch.

## Quick guide

1.  
Preparation

2.  
Place the sheet

3.  
Turn it to the left (90°)

4.  
Turn it to the center  
(180°)

5.  
Activate the cylinder

6.  
Heater on turn left

## Technical characteristics

Height: ..... 410 mm  
Width: ..... 545 mm  
Length: ..... 260 mm  
Weight: ..... 15 kg  
Power supply: ..... 220~240 V, 50/60 Hz  
Power consumption: ..... 450 W  
Pressure: ..... 4~6 bar