ACRYPOL FAST

RESINA ACRILICA TERMOPOLIMERIZZANTE RAPIDA.

FAST HEAT CURING ACRYLIC RESIN.

RESINA ACRÍLICA TERMOPOLIMERIZÁVEL RÁPIDA.

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(GB) FAST HEAT CURING ACRYLIC RESIN FOR COMPLETE, PARTIAL END COMBINED REMOVABLE DENTAL PROSTHESES AND RELINING.

DIRECTIONS: Fast heat curing acrylic resin liquid and powder with a high molecular weight for total and/or partial removable protheses and for combined protheses and relinings.

CONTRAINDICATIONS:

- Direct contact of non-polymerised material in the intraoral area.
- In case of proven allergy to components of Acry Pol Fast.
- All uses not indicated.

PREPARATION: Pour the desired quantity of liquid into a glass container, add the powder by tilting the container and tapping it on the support surface until a thin dry layer of powder remains on the surface. Wait for 20-30 seconds, remove any excess powder, mix for 30 seconds, cover and leave to rest. (Powder / liquid proportions 3: 1 by volume or 2: 1 by weight).
N.B.: At a temperature of 22°/23°C., Acry Pol Fast has a rest time of approximately 30 minutes. After the mixture is detached from the container using a spatula, the flask can be filled. The mixture retains its plasticity for approximately 15 minutes. The temperature has a considerable impact on rest time and resin plasticity, leaving the end result unchanged.

FLASK PREPARATION: After eliminating any traces of wax and washing thoroughly with boiling water, let the flask cool to a temperature between 20°C and a maximum of 40°C, and insulate the plaster with Acry Film or similar insulating material. Next, pack the prothesis mould with an adequate quantity of mixture. The two parts of the flask are then joined together with one or two pressing operations (at least 2 minutes), using a hydraulic press at a pressure of 120/150 kg/cm².

RECOMMENDED POLYMERISATION: Place the flask in boiling water, quickly bring the water back to the boil and simmer for 20 minutes. Remove the flask from the water and leave it to cool at room temperature. In the case of a high-thickness prothesis, immerse the flask in boiling water and turn off the device for 15 minutes. Then, restart the device and bring the water to the boil and keep the temperature constant for 20 minutes. Remove from the water and cool to room temperature.

FINISHING AND POLISHING: Finish as usual, with burs, rubber points, brushes and with pumice.

Mirror polish with the help of brushes and polishing products (Acry Lux Crema, Acry Lux Liquido, Acry Lux Pasta).

N.B.: Observing the above procedures will result in a denture with high stability and excellent tensile strength.

POWDER: Store in a cool and dry place

LIQUID: Store at a maximum temperature of 30°C

The product is designed for professional use in the dental sector. Follow the instructions for use. The manufacturer shall accept no responsibility for damages arising from different or inadequate use. The user must personally check the suitability of the product for the intended use, particularly if these uses are not reported in the instructions for use. Do not use after expiration.

Warning: Any serious accident that may occur with the device must be reported by the user to the manufacturer and the competent authority of the Member State where the user and/or patient are established.

LIQUID: Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid breathing vapours. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. Store in a well-ventilated place. Keep cool. Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains methyl methacrylate, ethylene dimethacrylate.

POWDER: Do not breathe dust. IF ON SKIN: Wash with plenty of soap and water.

Contains benzoyl peroxide, methyl methacrylate. May produce an allergic reaction.

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