ACRYC&B

ACRYC&B LIQUIDO **SELF CURING**
ACRYC&B LIQUIDO **REGULAR**
ACRYC&B LIQUIDO **SPECIAL**

ACRYC&B SMALTO
ACRYC&B COLLETTA
ACRYC&B DENTINA

RESINA ACRILICA PER RIVESTIMENTO ESTETICO DI PONTI E CORONE E PROVVISORI
ACRYLIC RESIN FOR AESTHETIC VENEERING OF CROWNS AND BRIDGES AND TEMPORARY FIXTURES
RESINA ACRÍLICA PARA REVESTIMENTO ESTÉTICO DE PONTES, COROAS E DENTES TEMPORÁRIOS
RÉSINE ACRYLIQUE POUR REVÊTEMENT ESTHÉTIQUE DE PONTS, COURONNES ET PROVISOIRES
RESINA ACRÍLICA PARA RECUBRIMIENTO ESTÉTICO DE PUENTES, CORONAS Y DIENTES PROVISIONALES
ACRYLHARZ ZUR ÄSTHETISCHEN BESCHICHTUNG VON BRÜCKEN UND KRONEN UND PROVISORIEN

(I) Istruzioni d'uso
(GB) Instructions for Use
(P) Instruções de uso
(F) Mode d'emploi
(E) Instrucciones de uso
(D) Gebrauchsanweisung
(CZ) Návod k použití
(LV) Lietošanas instrukcija
(TR) Kullanım için talimatlar
(SK) Inštrukcie na používanie
(RO) Instrucțiuni de folosire
(GR) Οδηγίες χρήσης
(BG) Инструкции за употреба
(RU) Инструкции по использованию

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(GB) ACRYLIC RESIN FOR AESTHETIC VENEERING OF CROWNS AND BRIDGES AND TEMPORARY FIXTURES

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- Repairs and applications of facing applications in temporary works.
- Fixing teeth for prostheses on frameworks.
- Production of temporary crowns and bridges.
- Application of acrylic facing applications (e.g. Acry Cette)

B) Acry C&B Liquido Regular (liquid) + Acry C&B SMALTO/COLLETTTO/DENTINA (powder) - Traditional muffle technique:

- Aesthetic covering for crowns and bridges.
- Production of temporary crowns and bridges.

C) Acry C&B Liquido Special (liquid) + Acry C&B SMALTO/COLLETTTO/DENTINA (powder) - Direct modelling technique with spatula or brush, with polymerisation in a special pressure device:

- Aesthetic veneering for crowns and bridges.
- Production of temporary crowns and bridges.

CONTRAINDICATIONS: - All unforeseen uses - Ascertained allergy to the components of Acry C&B - Direct use in the oral cavity - Direct contact of incorrectly polymerised material in the intraoral area.

PROCESSING:

A) Acry C&B Liquido Self Curing (liquid) + Acry C&B SMALTO/COLLETTTO/DENTINA (powder) - Cold working technique. Recommended technique for the construction of single temporary crowns or bridges, fixing teeth on frameworks, and to alter-characterise artificial resin teeth. Mix the powder with the Self-Curing Liquid, until a uniform mixture is obtained. The mixture is now ready for pouring; if you want a more plastic mix that is malleable with a spatula, let it rest for 1-2 minutes. Apply the material, filling the prepared templates or model with a spatula; if the material drips too much, add a little dry powder to the product to saturate the excess liquid. At a room temperature of 20°C the workability time is about 8 minutes and the mixture will harden after 15 - 18 minutes. We recommend polymerising in a pan or in a pressure polymeriser, with water at 45°/50°C for 10 minutes at a pressure of at least 2 atm.

FINISHING AND POLISHING: Finish as usual, with burs, rubber pads, brushes and pumice. Polish to a mirror shine using brushes and polishing products (Acry Lux Cream, Acry Lux Liquid, Acry Lux Paste) avoiding excessive overheating.

B) Acry C&B Liquido Regular (liquid) + Acry C&B SMALTO/COLLETTTO/DENTINA (powder) - Traditional muffle working technique. The traditional muffle working technique with hot counter-mould allows controlled stratification of the neck, dentin and enamel, with ample space for corrections. In this case, opacity (Ruthinium Opaker) the metal structure after having carefully degreased it.

MIXING Pour an adequate quantity of Regular Liquid into a container and saturate it with the powder. Mix vigorously with a spatula, cover the jar with a lid and let it rest for a few minutes. When it no longer remains attached to the spatula, the paste is ready for filling into a flask.

DENTINE STRATIFICATION: Apply the mixture starting from the incisal part under the stump if it is a jacket crown, distributing it with a spatula. Cover with a polyethylene sheet and close with the counter-mould, which has been heated in boiling water. Press slowly and let it rest for about three minutes. Reopen the flask, remove the polyethylene sheet and, with a blade, remove the excess dentin.

CERVICAL STRATIFICATION: With the same blade cut the dentin in the area of the collar, removing a small mass of material. Take care not to damage the opacifier if layering on a metal structure. Apply the collar mass, cover with a polyethylene sheet, and proceed with the second pressing, after heating the counter-mould. After reopening the flask, remove the excess material.

ENAMEL STRATIFICATION: With the same procedure used for the collar, remove the dentin mass in the incisal area, in greater quantities in the mesial and distal areas. Apply the enamel, cover with a polyethylene sheet, and proceed with the third pressing, with the same technique indicated in the previous steps. Before moving on to the final polymerisation phase, it is advisable to reopen the muffle to eliminate any burrs, apply a moistened polyethylene sheet and close definitively.

POLYMERISATION: Immerse the muffle in cold water; bring slowly to a boil in about 45 minutes and leave at this temperature for 30 minutes. Cool gradually and slowly.

FINISHING AND POLISHING: Finish as usual, with burs, rubber pads, brushes and pumice. Polish to a mirror shine using brushes and polishing products (Acry Lux Cream, Acry Lux Liquid, Acry Lux Paste) avoiding excessive overheating.

C) Acry C&B Liquido Special (liquid) + Acry C&B SMALTO/COLLETTTO/DENTINA (powder) - Direct modelling technique with spatula or brush, with polymerisation in a special pressure device. With this technique various traditional operating phases are eliminated, such as wax modelling, muffle setting and pressing, with a considerable saving of time.

MIXING: Once the opacifier is applied, if it is a prosthesis with a metal structure, or the insulator in the case of a jacket crown, pour an adequate amount of Special Liquid into a container and add resin, mixing until a good plastic mixture is obtained.

APPLICATION OF THE MASSES: After about 8 minutes, at a room temperature of 20°C, the mixture can be applied with a spatula or brush. The material remains workable for 15 - 20 minutes. For extensive work it is advisable to mix each single resin (collar, dentin, enamel) after completing the application of the previous one. First of all, apply the collar mass, blending it from the cervical part towards the centre of the tooth. Then apply the dentin, blending it towards the incisal area to leave room for the enamel mass.

POLYMERISATION: Put the product in the polymeriser set at 120°C and at least 6 atm for 15 minutes. If the plaster model also has to be put in the polymeriser, set the temperature to 100°C and pressure to 6 atm, for 30 minutes.

FINISHING AND POLISHING: Finish as usual, with burs, rubber pads, brushes and with pumice. Polish to a mirror shine using brushes and polishing products (Acry Lux Cream, Acry Lux Liquid, Acry Lux Paste) avoiding excessive overheating.

POWDER: Store in a dry place. Keep cool.

LIQUID: Store at temperatures not exceeding 30°C.

The product is designed for professional use in the dental sector. Follow the instructions for use. The manufacturer shall accept no responsibility for damages arising from different or inadequate use. The user must personally check the suitability of the product for the intended use, particularly if these uses are not reported in the instructions for use. Do not use after expiration. Warning: Any serious accident that may occur with the device must be reported by the user to the manufacturer and the competent authority of the Member State where the user and/or patient are established.

LIQUID: Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid breathing vapours. Wear protective gloves/protective clothing/eye protection/face protection. **IF ON SKIN:** Wash with plenty of soap and water. Store in a well-ventilated place. Keep cool. Dispose of contents/container in accordance with local/regional/national/international regulations. Contains methyl methacrylate, triethylene glycol dimethacrylate, ethylene dimethacrylate, N-N-dimethyl-p-toluidine.

POWDER: Contains benzoyl peroxide, methyl methacrylate. May produce an allergic reaction.

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Ruthinium® - Dental Manufacturing S.p.A.