## R-080030

E .....	ELECTROVIBRADORES .....	3
GB .....	ELECTRO-VIBRATORS .....	4
F .....	ELECTRO-VIBRATEURS .....	5
P .....	ELECTRO-VIBRADORES .....	6
I .....	ELETTROVIBRATORE .....	7

**MESTRA®**

TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60  
Tfno. + 34 94 453 03 88 Fax + 34 94 471 17 25  
E-mail: [mestra@mestra.es](mailto:mestra@mestra.es) - [www.mestra.es](http://www.mestra.es)  
48150 SONDIKA - BILBAO - ESPAÑA

Rev. 05/10/09

- |   |  |
|---|--|
| A | Tacos de goma / <b>Rubber bases</b> / <i>Ventouses</i> / Cunhas de borracha / <b>Tasselli di gomma</b>   |
| B | Regulador de vibraciones / <b>Vibration regulator</b> / <i>Régulateur de vibrations</i> / Regulador de vibrações / <b>Regolatore vibrazioni</b>              |
| C | Plataforma de vibración / <b>Vibration platform</b> / <i>Plateau de vibration</i> / Plataforma de vibração / <b>Piattaforma di vibrazione</b>                |
| D | Muesca apoyo cubetas / <b>Bowl support notch</b> / <i>Encoche d'appui pour porte empreintes</i> / Entalhadura apoio bacias / <b>Tacca appoggio vaschette</b> |

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The R-080030 is a practical electro-vibrator unit especially designed for compacting materials for the jewellery and dental mechanics sectors. To get the best performance from the unit we recommend that you read the following instructions carefully.

1. Place the unit on a plane, firm and perfectly horizontal surface. Do not place the unit near devices that could receive the vibrations. If you need reduce noise and vibrations, place the unit on a cloth or spongy material.
2. Connect the plug to a 220 V, 50/60 Hz, earthed power outlet.
3. Place the cylinder or elements to be vibrated on the unit's vibration platform. With bowls for impressions the mixture can also be compacted by resting the bowl firmly on the notch (D) in the vibration platform and holding down the bowl handle with your hand.
4. The on switch of the R-080030 is built into the regulator potentiometer (B). Turn the knob clockwise until you hear a click, and the electro-vibrator will start up.
5. Slowly move the unit's regulator (B) clockwise and check that the oscillation rate gradually increases. Set the regulator to a suitable amplitude for the work to be done.
6. After a short while, once the mixture is properly compacted, disconnect the unit from the mains by turning the regulator (B) anticlockwise until it clicks.

- ◆ Do not allow children or unskilled personnel to handle the unit.
- ◆ Do not pull, twist or bend the mains cable too far.
- ◆ Before cleaning the unit ensure that it is unplugged.
- ◆ Never immerse the unit in liquid when cleaning it.
- ◆ Do not place hot containers or objects on the vibration platform, as the rubber may be damaged.

#### TECHNICAL CHARACTERISTICS

Voltage:	220 V, 50/60 Hz
Work surface:	Ø 120 mm
Height:	100 mm
Weight:	2.5 kg
Power:	140 W
Vibrations per minute:	3000

#### MAINTENANCE AND CLEANING

To ensure a long, trouble-free lifetime and performance of the unit, we recommend as follows:

- After using the unit, always clean away any splashes of mixture before they dry.
- Clean the vibration platform and outer housing of the electro-vibrator regularly with a damp cloth. The rubber cover protecting the platform can be easily removed to facilitate cleaning.
- As a precaution, always disconnect the unit from the mains before cleaning it.

#### PRECAUTIONS

- ◆ Before connecting the unit ensure that the power supply is 220 V, 50/60 Hz and earthed.
- ◆ Do not use the unit in wet places or close to heat sources.
- ◆ Do not place easily broken containers or other objects on the vibration platform.

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