


# MESTRA®

BAÑO ULTRASONIDO "CALYPSO"  
"CALYPSO" ULTRASONIC CLEANER  
BAC À ULTRASONS "CALYPSO"  
R- 080422



**MESTRA®**  
**BAÑO ULTRASONIDO "CALYPSO"**  
**"CALYPSO" ULTRASONIC CLEANER**  
*BAC À ULTRASONS "CALYPSO"*  
R-080422

Nº	Referencia	DESCIPCION	DESCRIPTION
1	080420-14	Vaso de precipitados	Glass
2	080420-04	Junta tórica de vaso	Glass O-ring
4	080422-02	Tapa	Lid
5	080422-03	Cestillo	Gasket
6	080422-04	Porta vasos	Glass holder
7	080422-05	Timer	Timer
8	000250	Cable	Cable
9	080422-06	Carcasa	Housing
10	080422-07	Teclado	Keyboard
11	080422-08	Filtro de desagüe	Drain filter
12	080422-09	Cuba completa	Complete tank
13	080422-10	Chapas sujeta cuba	Tank holding sheets
14	080422-11	Tubo de desagüe	Drain pipe
15	080242-26	Contactador de cuba	Tank contacter
16	080422-13	Transductores	Transducers
17	080422-14	Muelle	Spring
18	080422-15	Control	Control board
19	080422-16	Resistencia	Resistor
20	080490-05	Fusible	Fuse
21	080030-21	Porta-fusible	Fuse holder
22	080422-17	Interruptor	Switch
23	080422-18	Filtro de red	Electric filter
24	080422-19	Codo	Elbow
25	080422-20	Llave de desagüe	Drain tap
26	080422-21	Racor salida	Outlet connector
27	080422-22	Chapa trasera	Back sheet
28	080422-23	Disipador	Heat sink
29	080070-22	Contactador disipador	Heat sing contacter
30	080490-31	Patas	Legs