

Alginate



Alginate for high precision dental impressions.
Available in both Chromatic and Dust Free. Packaging
of 24 bags of 450g each.

Chromatic:

Alginate with colour changes.

Flavour: mint.

Chromatic variations: violet – Mixing.
cyclamen – Working.
green – Mouth insertion.

Dust Free:

Fast mixing and setting time.

Flavour: mint.

Colour: green.

code

2040003	Alginate - Chromatic	450 g
2060003	Alginate - Dust Free	450 g

Chromatic Alginate

Mixing time	60"
Total Working Time	1'50"
Setting time	3'35"
First colour change	1'00"
Second colour change	1'40"
Reproduction of details	50μ
Recovery from deformation	95,5%
Strain in compression	12%
Compressive strenght	0,7Mpa
Flavour	Mint
Final colour	Green

Dust Free Alginate

Mixing time	60"
Total Working Time	1'50"
Setting time	3'35"
Reproduction of details	50μ
Recovery from deformation	95,5%
Strain in compression	12%
Compressive strenght	0,7Mpa
Flavour	Mint
Final colour	Green